



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

*Technically  
Speaking*

**CEDARLANE®**  
[www.cedarlanelabs.com](http://www.cedarlanelabs.com)



Conveniently Delivering You Today's Innovations  
for the Science of Tomorrow™

**Monoclonal Mouse antibody to Human Tyrosinase  
Pre-Diluted Antibody**

**CLSG90481**

**Lot:**

**Clone:** T311

**Immunogen:** BALB/C mice were injected with a recombinant tyrosinase protein

**Isotype:** IgG2a

**Size:** 6ml

**Format:** This antibody has been pretitered and quality controlled to work on formalin-fixed, paraffin embedded and acetone fixed cryosections tissue sections. No further titration is required.

**Staining Protocol:** Incubate for 30 minutes at room temperature. Optimal incubation conditions should be determined by the end user based upon the fixation conditions and staining system employed. High temperature treatment of formalin-fixed, paraffin embedded tissue sections using 1mM EDTA, pH 8.0 must be performed prior to the immunostaining.

**Specificity:** This antibody is specific to a cluster of proteins between 70-80 kD known as tyrosinase. This antibody does not show any cross-reaction with MAGE-1 and tyrosinase-related-protein-1, TRP-1/gp75. Staining of melanomas with this antibody showed tyrosinase in melanotic and amelanotic variants.

**Positive Control:** Melanoma

**Cellular Localization:** Cytoplasmic

**Storage:** Store at 2-8°C, do not use beyond expiration date stated on label.

**References:** Chen et al. Proc Natl Acad Sci USA 92:8125, 1995.

**Laboratory Reagent for Research Use Only**

RH/10/20/11

Visit our website for your local distributor.

**CEDARLANE®**



[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

An ISO 9001:2000 and ISO 13485:2003  
registered company.

**In CANADA: Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020  
e-mail: [general@cedarlanelabs.com](mailto:general@cedarlanelabs.com)

**In the USA: Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138  
e-mail: [service@cedarlanelabs.com](mailto:service@cedarlanelabs.com)