



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Recombinant Varicella Zoster Virus gE

CLVZV231
CLVZV231-2
CLVZV231-3

Introduction: VZV is closely related to the herpes simplex viruses(HSV), sharing much genome homology. The known envelope glycoproteins (gB, gC, gE, gH, gI, gK, gL) correspond with those in HSV, however there is no equivalent of HSV gD. VZV virions are spherical and 150-200 nm in diameter. Their lipid envelope encloses the nucleocapsid of 162 capsomeres arranged in a hexagonal form. Its DNA is a single, linear, double-stranded molecule, 125,000 nt long. The virus is very susceptible to disinfectants, notably sodium hypochlorite. Within the body it can be treated by a number of drugs and therapeutic agents including aciclovir, zoster-immune globulin (ZIG), and vidarabine.

Description: The E.Coli derived recombinant protein contains the VZV gE immunodominant regions, amino acids 48-135 and fused to a GST-Tag at N-terminus.

Purification Method: Varicella was purified by proprietary chromatographic technique.

Presentation: 100 µg (CLVZV231) or 0.5 mg (CLVZV231-2) or 1 mg (CLVZV231-3). 60mM NaCl, 50mM Tris-HCl pH 8.0, 0.25% Sarkosil, 10mM glutathione and 50% glycerol.

Stability: Varicella Protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

Purity: Varicella protein is >90% pure as determined by SDS-PAGE.

Specificity: Immunoreactive with sera of VZV-infected individuals.

Applications: ELISA, Western blots and Flow-through.

Laboratory Reagent For Research Use Only

MW 02/09/21

Visit our website for your local distributor.

CEDARLANE®



www.cedarlanelabs.com

An ISO 9001:2015 and ISO 13485:2016
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com