



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Catalog Number: 4442

Product Description: Corning® Elplasia™ Microcavity Plate, 96-Well, Black Clear Round Bottom, 500/400, Ultra-Low Attachment with Lid

**Component Materials:**

Plate walls - Virgin Polystyrene, black concentrate  
Plate bottom - Virgin Polystyrene  
Lid - Virgin Polystyrene  
Surface - Proprietary hydrogel

**Product Dimensions:**

Length of plate	-	5.01 in.	Working Volume	-	0.1-0.2mL
Width of plate	-	3.37 in.	Tolerances	-	+/- 0.05 in.
Height of plate with lid	-	0.68 in.			

**Microcavity Nominal Dimensions:**

Bottom diameter - 500µm  
Height - 400µm

**Sterilization:**

The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137 *Sterilization of healthcare products-Requirements for validation and routine control-Radiation sterilization*. Sterility Assurance Level: SAL 10<sup>-3</sup>

**Pyrogens:**

The product has been tested and has met the criteria established in the current version of ANSI/AAMI ST 72: *Bacterial Endotoxins - Test methodologies, routine testing, and alternative to batch testing*. Results: less than 0.1EU/mL.

**Tissue Culture:**

The product has been tested for the ability to promote the growth and formulation of 3D cellular spheroids of mammalian suspension cells in culture.

**Performance Testing:**

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Visual Attributes:	Visual examination of the product.
Packaging:	Inspection for seal and barrier integrity, accurate labeling, and correct product configuration.

**Lot Number Designation:**

8 Digit Lot Number: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

**Revision Date:**

03-17-20

Rev No: 2