



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com


[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

CORNING

Trace Elements Formulation

Cat. No.	25-021	25-022	25-023
Description	(A) 1000x	(B) 1000x	(C) 1000x
Units	mg/L	mg/L	mg/L

Components			
CuSO ₄ • 5H ₂ O	1.60	--	--
ZnSO ₄ • 7H ₂ O	863.00	--	--
Selenite • 2Na	17.30	--	--
Ferric citrate	1155.10	--	--
MnSO ₄ • H ₂ O	--	0.17	--
Na ₂ SiO ₃ • 9H ₂ O	--	140.00	--
Molybdic acid, ammonium salt	--	1.24	--
NH ₄ VO ₃	--	0.65	--
NiSO ₄ • 6H ₂ O	--	0.13	--
SnCl ₂ (anhydrous)	--	0.12	--
AlCl ₃ • 6H ₂ O	--	--	1.20
AgNO ₃	--	--	0.17
Ba(C ₂ H ₃ O ₂) ₂	--	--	2.55
KBr	--	--	0.12
CdCl ₂	--	--	2.28
CoCl ₂ • 6H ₂ O	--	--	2.38
CrCl ₃ (anhydrous)	--	--	0.32
NaF	--	--	4.20
GeO ₂	--	--	0.53
KI	--	--	0.17
RbCl	--	--	1.21
ZrOCl ₂ • 8H ₂ O	--	--	3.22



www.corning.com/lifesciences/solutions

At Corning, cells are in our culture. In our continuous efforts to improve efficiencies and develop new tools and technologies for life science researchers, we have scientists working in Corning R&D labs across the globe, doing what you do every day. From seeding starter cultures to expanding cells for assays, our technical experts understand your challenges and your increased need for more reliable cells and cellular material.

It is this expertise, plus a 160-year history of Corning innovation and manufacturing excellence, that puts us in a unique position to offer a beginning-to-end portfolio of high-quality, reliable cell culture consumables.

For additional product or technical information, please visit www.cellgro.com or call 1.800.235.5476.

Mediatech, Inc.
A Corning Subsidiary
9345 Discovery Boulevard
Manassas, VA 20109
t 800.235.5476
t 703.471.5955
f 703.467.9851
www.cellgro.com

CORNING | **FALCON** | **AXYGEN** | **GOSSELIN** | **PYREX**

Corning cellgro products are manufactured by Mediatech, Inc.
For a listing of Corning trademarks, visit us at www.cellgro.com/about-us/trademarks.html
All other trademarks in this document are the property of their respective owners.