



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Streptavidin(Peroxidase/HRP conjugated)

Catalog No. E-AB-1043

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Product Information

Purification Method: Affinity purification

Conjugation: HRP

Concentration: 0.5 mg/mL

Applications: IHC, WB, ELISA

Dilution: Immunohistochemistry in frozen section: 1:100-500

Immunohistochemistry in paraffin section: 1:100-200

Western Blotting: 1:3000-10000(Enhanced chemiluminescent detection)

Western Blotting: 1:500-3000(DAB detection)

Direct ELISA: 1:10000-100000(TMB detection)

Optimal working dilutions must be determined by end user. Dilute before when it's needed, and the residual should be discarded.

Buffer: PBS with 50% glycerol, PH7.4

Storage: Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

Focus on your research
Service for life science

Applications: Activ- Activation; Block- Blocking; Separation- Cell Separation ; Cell Sep-Neg- Cell Separation by Negative Selection; FA- Functional Assay; Neut- Neutralization; Stim- Stimulation; FCM- Flow Cytometry; ICFM: Intracellular Staining for Flow Cytometry; WB- Western Blotting; IHC- Immunohistochemistry; IF- Immunofluorescence; IP- Immunoprecipitation