



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Normal Quail IgY (egg-derived)

[exalpha.com/products/normal-quail-igy-egg-derived-2/QuIgY-100](https://www.exalpha.com/products/normal-quail-igy-egg-derived-2/QuIgY-100)

Catalog number: **QuIgY-100**

Product Type	Control antibodies
Units	100 mg
Host	Quail
Application	ELISA Western Blotting

Source

IgY was purified from the egg yolks of unimmunized quails.

Product

Liquid quail IgY fraction (opaque, light yellow in color).

Formulation: Phosphate buffered saline with 0.075% sodium azide

Purification Method: Egg yolk-derived purified IgY (non affinity purified)

Concentration: See vial for concentration.

Applications

Optimal working dilutions should be determined for your particular assay conditions.

Storage

The antibody may be stored at +4°C. For prolonged storage prepare appropriate aliquots and store at -20°C. Repeated thawing and freezing should be avoided. Working dilutions should be stored at +4°C, not refrozen, and preferably used the same day. An aliquot should be thawed slowly at ambient temperature and centrifuged to remove precipitate. This may result in a thin lipid layer floating on top and/or a visible precipitate. Removal of either will not affect the performance or the concentration of the product.

Shipping Conditions: Ship at ambient temperature.

Caution

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. It may contain hazardous ingredients. Please refer to the Safety Data Sheets (SDS) for additional

information and proper handling procedures. Dispose product remainders according to local regulations. This datasheet is as accurate as reasonably achievable, but Exalpha Biologicals accepts no liability for any inaccuracies or omissions in this information.

References

1. Gitoee, A., Sadeghi, G. & Karimi, A. 2017. Combination Effects of Organic and Inorganic Chromium on Production Performance, Reproductive Response, Immune Status, and Maternal Antibody Transmission in Breeder Quails Under Heat Stress. *Biol Trace Elem Res* <https://doi.org/10.1007/s12011-017-1205-x>

Safety Datasheet(s) for this product:

EA_Sodium Azide