



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Candida albicans antibody (FITC)

Cat. No. GTX04837-06

|             |                                            |
|-------------|--------------------------------------------|
| Host        | Rabbit                                     |
| Clonality   | Polyclonal                                 |
| Isotype     | IgG                                        |
| Application | ICC/IF, ELISA, Double immunodiffusion, IHC |
| Reactivity  | Candida albicans                           |

Package  
500 µl

## APPLICATION

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution     | Recommended dilution |
|------------------------|----------------------|
| ICC/IF                 | Assay dependent      |
| ELISA                  | Assay dependent      |
| Double immunodiffusion | Assay dependent      |
| IHC                    | Assay dependent      |

Not tested in other applications.

## Product Note

Recognizes numerous proteins in a soluble C. albicans extract (IEP). Has not been absorbed and does crossreact with other yeasts. Negative against human serum, urine and spinal fluid.

## PROPERTIES

|               |                                                                                                                                                                                                                                                |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                                         |
| Buffer        | PBS, 0.1% BSA                                                                                                                                                                                                                                  |
| Preservative  | 0.1% Sodium azide                                                                                                                                                                                                                              |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. Protect from light. |
| Concentration | 4-5 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                                     |
| Immunogen     | Candida albicans, Type A (ATCC #32354)_x000D_                                                                                                                                                                                                  |
| Purification  | Protein A chromatography purified IgG fraction covalently coupled with high purity isomer I of fluorescein isothiocyanate. Care is taken to ensure complete removal of any free fluorescein from the final product.                            |
| Conjugation   | Fluorescein isothiocyanate (FITC) <a href="#">Wavelength</a>                                                                                                                                                                                   |



For full product information, images and publications, please visit our [website](#).

Date 2024 / 12 / 12 Page 1 of 2

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).