



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# SARS-CoV-2 (COVID-19) Spike S2 (T716I, S982A, D1118H Mutant) (ECD) protein, His tag

**Cat No. GTX136023-pro**

<b>Application</b>	WB
<b>Species</b>	SARS Coronavirus 2


Package

100 µg

## PROPERTIES

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	batch dependent (Please refer to the vial label for the specific concentration.)
<b>Region/Sequence</b>	SARS-CoV-2 Spike of QHD43416.1 (686-1,213 a.a) with T716I, S982A and D1118H mutations, 6 stabilizing mutations (HexaPro mutations, including F817P, A892P, A899P, A942P, K986P, and V987P), and His tag at the C-terminus
<b>Expression System</b>	HEK293 cells
<b>Purity</b>	>95%
<b>Conjugation</b>	Unconjugated
<b>Note</b>	For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

## DATA IMAGES




### GTX136023-pro WB Image

SARS-CoV-2 (COVID-19) Spike S2 (T716I, S982A, D1118H Mutant) (ECD) protein, His tag (GTX136023-pro, 0.2 µg) was separated by 12% SDS-PAGE, and the membrane were blotted with SARS-CoV-2 (COVID-19) Spike S2 antibody [HL1038] (GTX635910, left) diluted at 1:20000, SARS-CoV-2 (COVID-19) Spike S2 antibody [HL1039] (GTX635911, middle) diluted at 1:20000, and SARS-CoV / SARS-CoV-2 (COVID-19) spike antibody [1A9] (GTX632604, right) diluted at 1:5000.



For full product information, images and publications, please visit our [website](#).



#### GTX136023-pro Image

1 $\mu$ g of GTX136023-pro SARS-CoV-2 (COVID-19) Spike S2 (T716I, S982A, D1118H Mutant) (ECD) protein, His tag was analyzed using SDS-PAGE and stained with coomassie blue and captured by monochrome camera.



For full product information, images and publications, please visit our [website](#).