



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

FluoroGel mounting medium

Cat No. GTX28214

Application ICC/IF, IHC-P, IHC-Fr, IHC

Reference (13)
 Package
 20 ml

PRODUCT

Product Description FluoroGel mounting medium - fluorescence

APPLICATION

Application Note

Prior to use turn the bottle dropper tip to the open position. Close the dropper tip tightly when finished mounting to avoid evaporation and solidifying of mounting medium.

PROPERTIES

Form Liquid

Storage Keep as concentrated solution, aliquot and store at 4°C.

Note For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).