



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# Immuno In Situ Mount™

**Cat No. GTX30924**

**Application** IHC

**Package**  
30 ml

## PRODUCT

**Product Description** Immuno In Situ Mount™

## APPLICATION

### Application Note

Apply 3-4 drops directly on the top of specimen and spread out evenly by tilting back and fourth or spread evenly with a 0.2 ml plastic pipette tip. Let stand at room temperature for about 60-90 minutes or heat slides at 37-40oC for 40-50 minutes.

## PROPERTIES

**Form** Liquid

**Storage** Room temperature

**Note** For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).