



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Endogenous peroxidase blocking kit

Cat No. GTX30967**Application** ICC/IF, IHC-P, IHC**Reference (3)**
Package
15 ml**PRODUCT****Product Description** Endogenous peroxidase blocking kit**APPLICATION****Application Note**

We recommend that endogenous peroxidase activity should be blocked following deparaffinization but prior to antigen retrieval.

PROPERTIES**Storage** store at 2-8 °C, DO NOT FREEZE**Note** For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.**DATA IMAGES****GTX30967 Image**

10µg of GTX109135-pro Galectin 2 recombinant protein analyzed using SDS-PAGE and stained with coomassie blue and captured by black and white camera.



For full product information, images and publications, please visit our [website](#).