



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com



www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Certificate of Analysis

LOT 6465

Product Information

Product: Pooled Normal Plasma
Product No.: 0010
Lot: 6465
Number of Donors: 56 Donors (28 Males 28 Females)
Expiry:  May 31, 2021
Recommended Storage:  -70° C. or colder

Product Description

George King Pooled Normal Plasma (PNP) is platelet poor fresh frozen citrated plasma pooled from human prescreened normal donors. No buffers or stabilizers added. Meets CLSI, CAP and CLOT-ED™ recommendations.

Donor selection criteria:
 Healthy donors
 Ages 18-65
 Non-Medicated
 Non-Smoking

Plasma collected by plasmapheresis in FDA licensed blood centers

Plasma tested non-reactive for HBsAg and HIV 1/2
 Anti-HCV Negative



George King Bio-Medical products meet the requirements of the European In Vitro Diagnostics Directive. **CE**

Test	Result	Instrument	Reagent	Normal Range
Prothrombin Time (PT)	13.5 s	Diagnostica Stago STA-R Evolution®	Diagnostica Stago STA® Neoplastine® CI Plus	12.8 - 15.2 s
Prothrombin Time (PT)	11.5 s	ACL TOP® 500	HemosIL® RecombiPlasTin	10.0 - 13.3 s
Prothrombin Time (PT)	10.4 s	Siemens Sysmex® CA-1500	Dade® Innovin®	9.7 - 11.7 s
Activated Partial Thromboplastin Time (aPTT)	30.7 s	Diagnostica Stago STA-R Evolution®	Diagnostica Stago STA® PTT A	27.3 - 36.1 s
Activated Partial Thromboplastin Time (aPTT)	29.4 s	ACL TOP® 500	HemosIL® SynthASil	25.5 - 35.5 s
Activated Partial Thromboplastin Time (aPTT)	28.9 s	Siemens Sysmex® CA-1500	Dade® Actin® FSL	23.6 - 31.9 s
Fibrinogen (Claus method)	>250 mg/dl	Diagnostica Stago STA-R Evolution®	Diagnostica Stago STA® Fibrinogen	170 - 400 mg/dl
Factors II, V, VII, and X	>86 %	Diagnostica Stago STA-R Evolution®	Diagnostica Stago STA® Neoplastine® CI Plus	50 - 150 %
Factors VIII, IX, XI, and XII	>80 %	Diagnostica Stago STA-R Evolution®	Diagnostica Stago STA® PTT A	50 - 150 %



George King Bio-Medical, Inc.
 Overland Park KS 66210 USA
 Phone +1 (913) 469-5464 | Toll Free (800) 255-5108
 Email plasma@kingbiomed.com
www.kingbiomed.com

