



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

8015171/8015172

Instant HEPES Buffer

Cat-No.	8015171	8015172
Size	10 packs/box	50 packs/box

Form	Crystalline powder, white to off-white granules
pH	7.7±0.2 at 25°C
Storage	Store at RT, protect from direct light
Shelf life	Powder: at least 36 months when stored at RT

Product Description

Each pack of Instant HEPES Buffer produces 1 L of 1 × buffer (50 mM Tris, 50 mM HEPES and 2 mM EDTA, 0.1% SDS). Instant HEPES Buffer is commonly used for protein analysis and biochemistry applications.

Preparation Notes

1. Add the entire contents of one pouch of Instant HEPES buffer into a glass beaker.
2. Add 500 ml of ultrapure water into the beaker and gently stir until completely dissolved.
3. Transfer the HEPES buffer solution into a measuring cylinder and add ultrapure water to 1 L final volume, yielding 1 × HEPES buffer, pH 7.7.
4. Transfer to a screw-cap glass bottle or similar container for storage. Store at RT.

