



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7


[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## WT-161

Cat. No.:	HY-100871
CAS No.:	1206731-57-8
Molecular Formula:	C <sub>27</sub> H <sub>30</sub> N <sub>4</sub> O <sub>3</sub>
Molecular Weight:	458.55
Storage:	Powder    -20°C    3 years 4°C    2 years In solvent   -80°C    2 years -20°C    1 year



### SOLVENT & SOLUBILITY

In Vitro

DMSO : ≥ 100 mg/mL (218.08 mM)

\* "≥" means soluble, but saturation unknown.

	<div>Solvent Concentration</div>	Mass	1 mg	5 mg	10 mg
Preparing Stock Solutions	1 mM		2.1808 mL	10.9039 mL	21.8079 mL
	5 mM		0.4362 mL	2.1808 mL	4.3616 mL
	10 mM		0.2181 mL	1.0904 mL	2.1808 mL

Please refer to the solubility information to select the appropriate solvent.

In Vivo

1. Add each solvent one by one: 10% DMSO >> 40% PEG300 >> 5% Tween-80 >> 45% saline  
Solubility: ≥ 2.5 mg/mL (5.45 mM); Clear solution

2. Add each solvent one by one: 10% DMSO >> 90% (20% SBE-β-CD in saline)  
Solubility: ≥ 2.5 mg/mL (5.45 mM); Clear solution

3. Add each solvent one by one: 10% DMSO >> 90% corn oil  
Solubility: ≥ 2.5 mg/mL (5.45 mM); Clear solution

### BIOLOGICAL ACTIVITY

Description	WT-161 is a potent and selective HDAC6 inhibitor with an IC <sub>50</sub> of 0.40 nM <sup>[1]</sup> . WT-161 also inhibits metallo-β-lactamase domain-containing protein 2 (MBLAC2) <sup>[2]</sup> .
IC <sub>50</sub> & Target	IC <sub>50</sub> : 0.40 nM (HDAC6), 51.61 nM (HDAC3) <sup>[1]</sup>
In Vitro	WT161 selectively inhibits HDAC6 and dramatically increases levels of acetylated α-tubulin (Ac-α-tubulin) with little effect on global lysine acetylation. WT161 induces significant toxicity in all multiple myeloma cell lines tested, with IC <sub>50</sub> s between 1.5 and 4.7 μM. WT161 in combination with bortezomib triggers significant accumulation of polyubiquitinated proteins and cell

stress, followed by caspase activation and apoptosis. More importantly, this combination treatment is effective in bortezomib-resistant cells and in the presence of bone marrow stromal cells, which have been shown to mediate multiple myeloma cell drug resistance<sup>[1]</sup>.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

#### In Vivo

WT161 shows toxicity at 100 mg/kg i.p., but WT161 is well tolerated at 50 mg/kg i.p.. Bortezomib combined with WT161 demonstrates a significant antitumor effect<sup>[1]</sup>.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

## PROTOCOL

#### Cell Assay <sup>[1]</sup>

MM.1S cells are treated with increasing concentrations of WT161 (0-10  $\mu$ M) for 48 hours. Cell viability is determined using the MTT assay<sup>[1]</sup>.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

#### Animal Administration <sup>[1]</sup>

Mice: Mice tumor xenograft are assigned into cohorts receiving vehicle (control), BTZ (0.5 mg/kg, i.v.), WT161 (50 mg/kg, i.p.), or BTZ+WT161. WT161 is administered for five consecutive days each week, and BTZ is given on a twice-weekly schedule. Caliper measurements of the longest perpendicular tumor diameters are performed on alternate days to estimate the tumor volume<sup>[1]</sup>.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

## CUSTOMER VALIDATION

- Cell Death Dis. 2023 Apr 6;14(4):250.
- J Mol Med (Berl). 2019 Aug;97(8):1183-1193.
- Biosci Rep. 2021 Apr 30;41(4):BSR20203905.
- J Appl Toxicol. 2023 Mar 1.
- Research Square Preprint. 2021 Jun.

See more customer validations on [www.MedChemExpress.com](http://www.MedChemExpress.com)

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: [tech@MedChemExpress.com](mailto:tech@MedChemExpress.com)

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA