

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



Proteins

Almitrine mesylate

Cat. No.: HY-107319 CAS No.: 29608-49-9 Molecular Formula: $C_{28}H_{37}F_{2}N_{7}O_{6}S_{2}$

Molecular Weight: 669.76

Target: Potassium Channel

Pathway: Membrane Transporter/Ion Channel

Storage: 4°C, sealed storage, away from moisture

* In solvent: -80°C, 6 months; -20°C, 1 month (sealed storage, away from moisture)

Product Data Sheet

SOLVENT & SOLUBILITY

In Vitro

DMSO: 100 mg/mL (149.31 mM; Need ultrasonic)

| Preparing Stock Solutions | Solvent Mass Concentration | 1 mg | 5 mg | 10 mg |
|------------------------------|-------------------------------|-----------|-----------|------------|
| | 1 mM | 1.4931 mL | 7.4654 mL | 14.9307 mL |
| | 5 mM | 0.2986 mL | 1.4931 mL | 2.9861 mL |
| | 10 mM | 0.1493 mL | 0.7465 mL | 1.4931 mL |

Please refer to the solubility information to select the appropriate solvent.

In Vivo

- 1. Add each solvent one by one: 10% DMSO >> 40% PEG300 >> 5% Tween-80 >> 45% saline Solubility: ≥ 2.5 mg/mL (3.73 mM); Clear solution
- 2. Add each solvent one by one: 10% DMSO >> 90% corn oil Solubility: ≥ 2.5 mg/mL (3.73 mM); Clear solution

BIOLOGICAL ACTIVITY

| Description | Almitrine mesylate, a peripheral chemoreceptor agonist, inhibits selectively the Ca ²⁺ -dependent K ⁺ channel. | |
|---------------------------|--|--|
| IC ₅₀ & Target | $K^+channel^{[1]}$ | |
| In Vitro | Almitrine inhibits the activity of a high-conductance (152±13 pS), Ca^{2+} -dependent K ⁺ channel by decreasing its open probability. The IC ₅₀ value of the effect is 0.22 μ M. The inhibitory effect of Almitrine on Ca^{2+} -dependent K ⁺ channels also is observed in GH3 cells. Almitrine at concentrations up to 10 μ M does not affect whole-cell voltage-dependent K ⁺ , Ca^{2+} , or Na ⁺ currents in rat or rabbit cells. However, this concentration of Almitrine significantly inhibits the Ca^{2+} -dependent component of K ⁺ currents in rat chemoreceptor cells ^[1] . MCE has not independently confirmed the accuracy of these methods. They are for reference only. | |

In Vivo

Almitrine acts via the peripheral arterial chemoreceptors raising carotid sinus nerve output and minute ventilation.

Almitrine also has a pulmonary vascular action causing a dose-dependent constriction and dilatation. At low doses Almitrine enhances hypoxic pulmonary vasoconstriction and may improve the overall ventilation/perfusion ratio^[2].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

PROTOCOL

Animal
Administration [2]

Rats^[2]

8-week-old male spf Wistar rats are used. Rats are anaesthetized with Thiobarbiturate inactin (BYK, 100 mg/kg, i.p.). The interaction of the ventilatory response to hypoxia and an intermittent (2 min on, 1 min break) low-dose (10 μ g/kg per min) and high-dose (80 μ g/kg per min) infusion of S9581 or Almitrine is tested in control and chronically hypoxic rats. S9581 or Almitrine is infused intravenously (100 μ g/ mL). Inspired oxygen levels were controlled by passing oxygen-nitrogen mixtures across the tracheal port at a flow rate of 3-41 min^{-1[2]}.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

REFERENCES

[1]. López-López JR, et al. Effects of Almitrine bismesylate on the ionic currents of chemoreceptor cells from the carotid body. Mol Pharmacol. 1998 Feb;53(2):330-9.

[2]. Bee D, et al. An analysis of the action of an analogue of Almitrine bismesylate in the rat model of hypoxic lung disease. Exp Physiol. 1992 Nov;77(6):819-28.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA