

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



Product Data Sheet

Fexagratinib

Cat. No.: HY-13330 CAS No.: 1035270-39-3 Molecular Formula: $C_{26}H_{33}N_5O_3$ Molecular Weight: 463.57 Target: FGFR

Pathway: Protein Tyrosine Kinase/RTK

Storage: Powder -20°C 3 years

4°C 2 years

In solvent -80°C 2 years

-20°C 1 year

SOLVENT & SOLUBILITY

In Vitro

DMSO: 125 mg/mL (269.65 mM; Need ultrasonic)

Preparing Stock Solutions	Solvent Mass Concentration	1 mg	5 mg	10 mg
	1 mM	2.1572 mL	10.7859 mL	21.5717 mL
	5 mM	0.4314 mL	2.1572 mL	4.3143 mL
	10 mM	0.2157 mL	1.0786 mL	2.1572 mL

Please refer to the solubility information to select the appropriate solvent.

In Vivo

- 1. Add each solvent one by one: 1% CMC-Na/saline water Solubility: 3.33 mg/mL (7.18 mM); Suspended solution; Need ultrasonic and warming and heat to 60°C
- 2. Add each solvent one by one: 10% DMSO >> 40% PEG300 >> 5% Tween-80 >> 45% saline Solubility: ≥ 2.5 mg/mL (5.39 mM); Clear solution
- 3. Add each solvent one by one: 10% DMSO >> 90% (20% SBE-β-CD in saline) Solubility: ≥ 2.08 mg/mL (4.49 mM); Clear solution
- 4. Add each solvent one by one: 10% DMSO >> 90% corn oil Solubility: ≥ 2.08 mg/mL (4.49 mM); Clear solution

BIOLOGICAL ACTIVITY

Description Fexagratinib (AZD4547; ADSK091) is a potent inhibitor of the FGFR family with IC₅₀s of 0.2 nM, 2.5 nM, 1.8 nM, and 165 nM for FGFR1, FGFR2, FGFR3, and FGFR4, respectively.

Page 1 of 3

In Vitro

Fexagratinib also inhibits recombinant VEGFR2 (KDR) kinase activity with an IC $_{50}$ of 24 nM. In KG1a, Sum52-PE, MCF7, and KMS11 cell lines, Fexagratinib potently inhibits autophosphorylation of FGFR1, 2, and 3 tyrosine kinases (IC $_{50}$ values of 12, 2, and 40 nM, respectively) and displays weaker inhibition of FGFR4 cellular kinase activity (IC $_{50}$ =142 nM). Significantly weaker inhibitory activity is observed versus cellular KDR and IGFR ligand-induced phosphorylation (IC $_{50}$ values of 258 and 828 nM, respectively), representing approximately 20- and 70-fold selectivity over cellular FGFR1. Besides, Fexagratinib potently inhibits FGFR phosphorylation and downstream signaling affected through FRS2, PLCy, and MAPK at the cellular level^[1]. MCE has not independently confirmed the accuracy of these methods. They are for reference only.

In Vivo

Female SCID mice bearing KMS11 tumors are randomized and treated chronically with Fexagratinib at a range of well-tolerated doses. Oral Fexagratinib treatment results in dose-dependent tumor growth inhibition. Twice daily administration of Fexagratinib at 3 mg/kg gives statistically significant tumor growth inhibition of 53% (P<0.0005 by one-tailed t test) when compare with vehicle-treated controls, whereas doses of 12.5 mg/kg once daily and 6.25 mg/kg twice daily results in complete tumor stasis (P<0.0001). A further efficacy study in the KG1a model with 12.5 mg/kg once daily Fexagratinib results in 65% tumor growth inhibition (P=0.002)^[1].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

PROTOCOL

Cell Assay [1]

Cell lines are incubated with fixed concentrations of AZD4547 for 72 hours. For fluorescence-activated cell sorting (FACS), cells are fixed with 70% ethanol and then incubated with propidium iodide/RNase A labeling solution. Cell-cycle profiles are assessed with a FACSCalibur instrument and CellQuest analysis software. For apoptotic analysis, cells and media are gently harvested and centrifuged, followed by washing of cell pellets. Cells are then processed for FITC staining and propidium iodide uptake. The proportion of cells staining positive for Annexin V are then assessed with a FACSCalibur instrument and quadrant sorting is done by CellQuest analysis software^[1].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

Animal Administration [1]

Mice^[1]

Swiss derived nude (nu/nu) and severe combined immunodeficient mice (SCID) are used. Tumor xenografts are established by s.c. injection into the left flank with 0.1 mL tumor cells (1×10^6 for LoVo, 1×10^7 for HCT-15, and 1×10^7 for Calu-6) or 0.2 mL (2×10^7 for KMS11 and KG1a) mixed 1:1 with Matrigel, with the exception of LoVo and HCT-15, which do not include Matrigel. Mice are randomized into control and treatment groups (AZD4547, 1.5-50 mg/kg, once daily or twice daily by oral gavage) when tumors reach the determined size of more than 0.2 cm³. Tumor volume (measured by caliper), animal body weight, and tumor condition are recorded twice weekly for the duration of the study. Tumor volume is calculated.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

CUSTOMER VALIDATION

- Bioact Mater. 2 January 2022.
- Sci Transl Med. 2018 Jul 18;10(450):eaaq1093.
- Nat Commun. 2020 Aug 13;11(1):4053.
- J Exp Clin Cancer Res. 2023 Apr 3;42(1):79.
- J Nanobiotechnology. 2023 Feb 17;21(1):55.

See more customer validations on www.MedChemExpress.com

REFERENCES

1]. Gavine PR, et al. AZD4547: ar 2045-2056.	n orally bioavailable, potent	, and selective inhibitor of the fi	broblast growth factor receptor tyro	sine kinase family. Cancer Res, 2012, 72(8),
Caution: Product has not been fully validated for medical applications. For research use only.				
	Tel: 609-228-6898 Address: 1	Fax: 609-228-5909 Deer Park Dr, Suite Q, Monm	E-mail: tech@MedChemEx outh Junction, NJ 08852, USA	press.com

Page 3 of 3 www.MedChemExpress.com