



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Cross-linked dextran G 150

Cat. No.: HY-156045B
CAS No.: 12774-36-6
Target: Biochemical Assay Reagents
Pathway: Others
Storage: Store at room temperature 3 years

Cross-linked dextran G 150

BIOLOGICAL ACTIVITY

Description

Cross-linked dextran G 150 is a hydrophilic gel that can be used as a gel filter filler (Particle size range: 40-120 μm ; Sphere protein separation range: 5k-300k Da)^[1].

REFERENCES

[1]. Ogston AG, et al. The thermodynamics of interaction between Sephadex and penetrating solutes. Biochem J. 1970 Jan;116(2):171-5.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA