



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

MEGDSRGG-(Ala-¹³C₃,¹⁵N)-GGEA-(Leu-¹³C₆,¹⁵N)-DPESAASSLSGAS TFA

Cat. No.:	HY-158994S
Molecular Formula:	C ₈₇ ¹³ C ₉ H ₁₅₆ N ₂₈ ¹⁵ N ₂ O ₄₄ S.xC ₂ HF ₃ O ₂
Sequence:	Met-Glu-Gly-Asp-Ser-Arg-Gly-Gly-(Ala-13C ₃ ,15N)-Gly-Gly-Glu-Ala-(Leu-13C ₆ ,15N)-Asp-Pro-Glu-Ser-Ala-Ala-Ser-Ser-Leu-Ser-Gly-Ala-Ser Met-Glu-Gly-Asp-Ser-Arg-Gly-Gly-(Ala-13C ₃ ,15N)-Gly-Gly-Glu-Ala-(Leu-13C ₆ ,15N)-Asp-Pro-Glu-Ser-Ala-Ala-Ser-Ser-Leu-Ser-Gly-Ala-Ser (TFA salt)
Sequence Shortening:	MEGDSRGG-{Ala-13C ₃ ,15N}-GGEA-{Leu-13C ₆ ,15N}-DPESAASSLSGAS
Target:	Isotope-Labeled Compounds
Pathway:	Others
Storage:	Please store the product under the recommended conditions in the Certificate of Analysis.

REFERENCES

[1]. Russak EM, et al. Impact of Deuterium Substitution on the Pharmacokinetics of Pharmaceuticals. Ann Pharmacother. 2019 Feb;53(2):211-216.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA