



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com


www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)



Pyrene-1,3,6,8-tetrasulfonic acid tetrasodium

Cat. No.:	HY-D0142
CAS No.:	59572-10-0
Molecular Formula:	$C_{16}H_6Na_4O_{12}S_4$
Molecular Weight:	610.43
Target:	Fluorescent Dye
Pathway:	Others
Storage:	4°C, sealed storage, away from moisture * In solvent : -80°C, 6 months; -20°C, 1 month (sealed storage, away from moisture)



BIOLOGICAL ACTIVITY

Description

Pyrene-1,3,6,8-tetrasulfonic acid tetrasodium is a fluorescent dye and pH indicator, also as a ligand of multifunctional metal-organic framework. Pyrene-1,3,6,8-tetrasulfonic acid tetrasodium has been used to detect CO₂ release^{[1][2][3]}.

REFERENCES

- [1]. Levin, P.P., et al. Kinetics of hydrated electron reactions with phosphate anions: a laser photolysis study. Russ Chem Bull 56, 1325–1328 (2007).
- [2]. Amalia García-García, et al. Slow relaxation of magnetization and luminescence properties of a novel dysprosium and pyrene-1,3,6,8-tetrasulfonate based MOF. New Journal of Chemistry. Issue 2, 2018.
- [3]. Bin Yan, et al. Manipulation of block copolymer vesicles using CO₂: dissociation or “breathing”. Soft Matter. Issue 6, 2013.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA