



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7


mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Direct Blue 86

Cat. No.:	HY-D0953
CAS No.:	1330-38-7
Molecular Formula:	C ₃₂ H ₁₄ CuN ₈ O ₆ S ₂ ·2Na
Molecular Weight:	780.16
Target:	Fluorescent Dye
Pathway:	Others
Storage:	RT, stored under nitrogen, away from moisture



BIOLOGICAL ACTIVITY

Description

Direct Blue 86 (Solvent Blue 38) is a myelin-sheath stain, commonly utilized in microscopy to detect demyelination in the central nervous system. Direct Blue 86 also is a dye with various applications including as a commercial dye in the printing of cotton and mucilage glue fabrics^[1].

REFERENCES

[1]. Ahmed El Nemr, et al. Removal of direct blue-86 from aqueous solution by new activated carbon developed from orange peel. J Hazard Mater. 2009 Jan 15;161(1):102-10.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA