



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Warm-Start M-MLV Reverse Transcriptase

Cat. No.:	HY-E70403	
Target:	Biochemical Assay Reagents	
Pathway:	Others	
Storage:	Please store the product under the recommended conditions in the Certificate of Analysis.	Warm-Start M-MLV Reverse Transcriptase

BIOLOGICAL ACTIVITY

Description

Warm-Start M-MLV Reverse Transcriptase is a reverse transcriptase active below 37°C and can completely release reverse transcription activity at 55°C without affecting cDNA synthesis^[1].

REFERENCES

[1]. Igor P. Osorbin, et al. M-MuLV reverse transcriptase: Selected properties and improved mutants. Comput Struct Biotechnol J. 2021; 19: 6315–6327.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA