



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Phalloidin

Cat. No.: HY-P0028

CAS No.: 17466-45-4

Molecular Formula: $C_{35}H_{48}N_8O_{11}S$

Molecular Weight: 788.87

A-(d-Thr)-C-(Hyp)-AWL (Sulfide bridge: Cys₃-Trp₆)

Storage: Sealed storage, away from moisture and light

Powder -80°C 2 years
 -20°C 1 year

* In solvent : -80°C, 6 months; -20°C, 1 month (sealed storage, away from moisture and light)

SOLVENT & SOLUBILITY

In Vitro

DMSO : 10 mg/mL (12.68 mM; Need ultrasonic)

	Solvent Concentration	Mass	1 mg	5 mg	10 mg
Preparing Stock Solutions	1 mM		1.2676 mL	6.3382 mL	12.6764 mL
	5 mM		0.2535 mL	1.2676 mL	2.5353 mL
	10 mM		0.1268 mL	0.6338 mL	1.2676 mL

Please refer to the solubility information to select the appropriate solvent.

BIOLOGICAL ACTIVITY

Description

Phalloidin is a mushroom-derived toxin which can be used to label F-actin of the cytoskeleton with fluorochrome^[1].

In Vitro

Phalloidin staining will not work with methanol-fixed cells, probably because of the disruption of actin filament integrity, but it will work with cells fixed in 0.2% glutaraldehyde in PBS^[1].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

PROTOCOL

Cell Assay ^[1]

Prepare a 1:200 dilution of the stock solution (300 units/mL) of fluorescein Phalloidin or rhodamine Phalloidin using PBS. Aspirate the cell medium from the cells grown on glass coverslips and rinse the cells three times with PBS. Fix the cells for 10 min in the 3.7% formaldehyde solution. Rinse the fixed cells three times for 5 min each in PBS. Permeabilize the cells for 5 min in the 0.2% Triton X-100 solution. Rinse the fixed cells three times with PBS. Label the cells with the fluorescein or rhodamine Phalloidin for 5 to 10 min at room temperature. Rinse the cells three times in PBS for 5 min each time. Mount the cell and view the cells using either the fluorescein or rhodamine filter set, depending on the probe^[1].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

CUSTOMER VALIDATION

- Drug Resist Updat. 2024 Mar, 73, 101056.
- Cell Commun Signal. 2024 Jan 17;22(1):44.
- J Mol Cell Biol. 2022 Jan 11;mjac001.
- J Cell Mol Med. 2024 Feb;28(4):e18130.
- Poult Sci. 23 February 2022, 101801.

See more customer validations on www.MedChemExpress.com

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA