



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

GelMA, 90% methacrylation

Cat. No.:	HY-W763582B	
CAS No.:	2280018-77-9	
Target:	Biochemical Assay Reagents	
Pathway:	Others	GelMA, 90% methacrylation
Storage:	Please store the product under the recommended conditions in the Certificate of Analysis.	

BIOLOGICAL ACTIVITY

Description

Gelatin Methacryloyl (GelMA) is a derivative obtained by the reaction of chitin anhydride (methacrylic anhydride, MA) and gelatin. The hydrogel formed by Gelatin Methacryloyl has good biocompatibility and biodegradability, photocrosslinking ability and adjustable physical properties^[1]. GelMA, 90% methacrylation needs to self-assemble into fibrous hydrogel under the action of photoinitiator LAP (HY-44076), and target bioactive adhesion sites, play an inherent supporting role for tissue cells and biodegradable activity.

Application: cell culture, biological 3D printing, tissue engineering, etc.

REFERENCES

[1]. Piao Y, et al., Biomedical applications of gelatin methacryloyl hydrogels[J]. Engineered Regeneration, 2021, 2: 47-56.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA