



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Sodium carboxymethyl cellulose (Viscosity:5000-15000 mPa.s)

Cat. No.:	HY-Y1889			Sodium carboxymethyl cellulose
CAS No.:	9004-32-4			
Target:	Biochemical Assay Reagents			
Pathway:	Others			
Storage:	Powder	-20°C 3 years 4°C 2 years		

SOLVENT & SOLUBILITY

In Vitro	H ₂ O : 6.25 mg/mL (Need ultrasonic)
-----------------	---

BIOLOGICAL ACTIVITY

Description	Sodium carboxymethyl cellulose (Viscosity:5000-15000 mPa.s) is the sodium salt of cellulose arboxymethyl and frequently used as viscous agent, paste and barrier agent.
--------------------	---

In Vivo	<p>Protocol for preparing 0.5% CMC-Na Solution Measure 0.5g of dry CMC-Na and dissolved in 100 ml ddH₂O/0.9% Saline (0.9 g NaCl in 100 ml ddH₂O) to make a clear solution. Under the condition of stirring and heating (50-65°C), adding CMC-Na slowly to ddH₂O/0.9% Saline helps to accelerate dissolution.</p> <p>Note</p> <ol style="list-style-type: none"> You must ensure that your CMC-Na solution does not exist solid-liquid separation phenomenon. The solution is in a uniform and transparent state has no particles in it. Completely dissolution of CMC-Na may requires 4 hours or more longer. <p>In a pharmacological test, CMC-Na (oral;5% in water; 1 year) is well tolerated in rats^[2]. In an acute oral toxicity study in female mice, LD₅₀ of CMC-Na for female mice is 14 g/kg body weight of mice, equivalent to 9.8 g/kg body weight of rat, categorized as practically non-toxic according to Loomis criteria (LD₅₀ 5-15g/kg body weight of rat)^[3]. MCE has not independently confirmed the accuracy of these methods. They are for reference only.</p>
----------------	---

REFERENCES

- [1]. Du XH, et al. Dexamethasone and sodium carboxymethyl cellulose prevent postoperative intraperitoneal adhesions in rats. *Braz J Med Biol Res.* 2015 Apr;48(4):344-8.
- [2]. Common Vehicles for Nonclinical Evaluation of Therapeutic Agents.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA