

# Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



# Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

## Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

## SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



## Product datasheet

# Mouse anti Basal cell Cytokeratin

(%) nordicmubio.com/products/mouse-anti-basal-cell-cytokeratin/MUB0312P-CE slash IVD

Catalog number: MUB0312P-CE/IVD

Clone	RCK103
Isotype	lgG1
Product Type	Primary Antibodies
Units	0.1 mg
Host	Mouse
Species Reactivity	Canine Chicken Guinea Pig Hamster Human Quail Rabbit Rat Swine Zebrafish
Application	Flow Cytometry Immunocytochemistry Immunohistochemistry (frozen)

Western Blotting

## **Background**

Cytokeratins are a subfamily of intermediate filament proteins and are characterized by a remarkable biochemical diversity, represented in Human epithelial tissues by at least 20 different polypeptides. They range in molecular weight between 40 kDa and 68 kDa and isoelectric pH between 4.9-7.8. The individual Human Cytokeratins are numbered 1 to 20. The various epithelia in the Human body usually express Cytokeratins which are not only characteristic of the type of epithelium, but also related to the degree of matuRation or differentiation within an epithelium. Cytokeratin subtype expression patterns are used

to an increasing extent in the distinction of different types of epithelial malignancies. The Cytokeratin antibodies are not only of assistance in the differential diagnosis of tumors using immunohistochemistry on tissue sections, but are also a useful tool in cytopathology and flow cytometric assays.

#### Source

RCK103 is a Mouse monoclonal IgG1 antibody derived by fusion of SP2/0-Ag14 Mouse myeloma cells with spleen cells from a BALB/c Mouse immunized with a mix of cell preparations containing Human Cytokeratins.

#### **Product**

Each vial contains 100 ul 1 mg/ml purified monoclonal antibody in PBS containing 0.09% sodium azide.

Formulation: Each vial contains 100 ul 1 mg/ml purified monoclonal antibody in PBS containing 0.09% sodium azide.

## **Specificity**

RCK103 is a Cytokeratin antibody recognizing (amongst others) Cytokeratin 5. This monoclonal antibody stains basal cells in combined and stRatified epithelial tissues. It recognizes the stem cell population, including the so-called amplifying cells in the prostate epithelium.

## **Applications**

RCK103 is useful for immunoblotting, immunocytochemistry, immunohistochemistry on frozen tissues and flow cytometry. Optimal antibody dilution should be determined by titration; recommended range is 1:100 – 1:200 for flow cytometry, and for immunohistochemistry with avidin-biotinylated Horseradish peroxidase complex (ABC) as detection reagent, and 1:100 – 1:1000 for immunoblotting applications.

### Storage

The antibody is shipped at ambient temperature and may be stored at +4°C. For prolonged storage prepare appropriate aliquots and store at or below -20°C. Prior to use, an aliquot is thawed slowly in the dark at ambient temperature, spun down again and used to prepare working dilutions by adding sterile phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided. Working dilutions should be stored at +4°C, not refrozen, and preferably used the same day. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance or the concentration of the product.

## Caution

When used for in vitro diagnostic purposes results must be put within the context of other diagnostic tests as well as the clinical history of the patient by a certified professional before final interpretation. Analyses performed with this antibody should be paralleled by positive and negative controls. If unexpected results are obtained which cannot be attributed to differences in laboratory procedures, please contact us. This product may

contain hazardous ingredients. Please refer to the Safety Data Sheets (SDS) for additional information and proper handling procedures. Dispose product remainders according to local regulations. This datasheet is as accurate as reasonably achievable, but Exalpha Biologicals accepts no liability for any inaccuracies or omissions in this information.

#### References

1. Feitz, W. F., Debruyne, F. M., Vooijs, G. P., Herman, C. J., and Ramaekers, F. C. (1986). Intermediate filament proteins as tissue specific markers in normal and malignant urological tissues, J Urol 136, 922-31. 2. Verhagen, A. P., Aalders, T. W., Ramaekers, F. C., Debruyne, F. M., and Schalken, J. A. (1988). Differential expression of Keratins in the basal and luminal compartments of Rat prostatic epithelium during degeneRation and regeneration, Prostate 13, 25-38. 3. Schaafsma, H. E., Ramaekers, F. C., van Muijen, G. N., Ooms, E. C., and Ruiter, D. J. (1989). Distribution of Cytokeratin polypeptides in epithelia of the adult Human urinary tract, Histochemistry 91, 151-9. 4. Smedts, F., Ramaekers, F., Robben, H., Pruszczynski, M., van Muijen, G., Lane, B., Leigh, I., and Vooijs, P. (1990). Changing patterns of Keratin expression during progression of cervical intraepithelial neoplasia, Am J Pathol 136, 657-68. 5. Kuijpers, W., Tonnaer, E. L., Peters, T. A., and Ramaekers, F. C. (1991). Expression of intermediate filament proteins in the mature inner ear of the Rat and Guinea Pig, Hear Res 52, 133-46. 6. Vos, J. H., van den Ingh, T. S., de Neijs, M., van Mil, F. N., Ivanyi, D., and Ramaekers, F. C. (1992). Immunohistochemistry with Keratin monoclonal antibodies in canine tissues: urogenital tract, respiratory tract, (neuro-) endocrine tissues, choroid plexus and spinal cord, J Vet Med 39, 721-40. 7. Vos, J. H., van den Ingh, T. S., Ramaekers, F. C., Molenbeek, R. F., de Neijs, M., van Mil, F. N., and Ivanyi, D. (1993). The expression of Keratins, vimentin, neurofilament proteins, smooth muscle actin, neuron-specific enolase, and synaptophysin in tumors of the specific glands in the canine anal region, Vet Pathol 30, 352-61. 8. Verhagen, A. P., Ramaekers, F. C., Aalders, T. W., Schaafsma, H. E., Debruyne, F. M., and Schalken, J. A. (1992). Colocalization of basal and luminal cell-type Cytokeratins in Human prostate cancer, Cancer Res 52, 6182-7. 9. van Leenders, G., Dijkman, H., Hulsbergen-van de Kaa, C., Ruiter, D., and Schalken, J. (2000). Demonstration of intermediate cells during Human prostate epithelial differentiation in situ and in vitro using triple-staining confocal scanning microscopy, Lab Invest 80, 1251-8.

## **CE Mark**

CE

### Protein Reference(s)

Database Name: UniProt

Accession Number: multiple

## Safety Datasheet(s) for this product:

NM Sodium Azide

