



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



---

## Mouse anti Cytokeratin 7 / Keratin K7

---

 [nordicmubio.com/products/mouse-anti-cytokeratin-7-keratin-k7/MUB0315P-CE\\_slash\\_IVD](https://nordicmubio.com/products/mouse-anti-cytokeratin-7-keratin-k7/MUB0315P-CE_slash_IVD)

Catalog number: **MUB0315P-CE/IVD**

Clone	RCK105
Isotype	IgG1
Product Type	Primary Antibodies
Units	0.1 mg
Host	Mouse
Species Reactivity	Canine Caprine Feline Hamster Human Mouse Rat Zebrafish
Application	Flow Cytometry Immunocytochemistry Immunohistochemistry (frozen) Western Blotting

### Background

Cytokeratins are a subfamily of intermediate filament proteins and are characterized by a remarkable biochemical diversity, represented in Human epithelial tissues by at least 20 different polypeptides. They range in molecular weight between 40 kDa and 68 kDa and isoelectric pH between 4.9 – 7.8. The individual Human Cytokeratins are numbered 1 to 20. The various epithelia in the Human body usually express Cytokeratins which are not only characteristic of the type of epithelium, but also related to the degree of maturation or differentiation within an epithelium. Cytokeratin subtype expression patterns are used in the distinction of different types of epithelial malignancies. The Cytokeratin antibodies

are not only of assistance in the differential diagnosis of tumors using immunohistochemistry on tissue sections, but are also a useful tool in cytopathology and flow cytometric assays.

### **Source**

RCK105 is a mouse monoclonal IgG1 antibody derived by fusion of SP2/0-Ag14 mouse myeloma cells with spleen cells from a BALB/c mouse immunized with cytokeratins from the human bladder carcinoma cell line T24.

### **Product**

Each vial contains 100 ul 1 mg/ml purified monoclonal antibody in PBS containing 0.09% sodium azide.

*Formulation:* Each vial contains 100 ul 1 mg/ml purified monoclonal antibody in PBS containing 0.09% sodium azide.

### **Specificity**

RCK105 reacts exclusively with Cytokeratin 7 which is present in a subgroup of glandular epithelia and their tumors, as well as transitional epithelium and transitional carcinoma.

### **Applications**

RCK105 is suitable for immunoblotting, immunocytochemistry on permeabilised cells, immunohistochemistry on frozen tissues and flow cytometry. Optimal antibody dilution should be determined by titration; recommended range is 1:100 – 1:200 for flow cytometry and immunohistochemistry with avidin-biotinylated Horseradish peroxidase complex (ABC) as detection reagent, and 1:100 – 1:1000 for immunoblotting applications.

### **Storage**

The antibody is shipped at ambient temperature and may be stored at +4°C. For prolonged storage prepare appropriate aliquots and store at or below -20°C. Prior to use, an aliquot is thawed slowly in the dark at ambient temperature, spun down again and used to prepare working dilutions by adding sterile phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided. Working dilutions should be stored at +4°C, not refrozen, and preferably used the same day. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance or the concentration of the product.

### **Caution**

When used for in vitro diagnostic purposes results must be put within the context of other diagnostic tests as well as the clinical history of the patient by a certified professional before final interpretation. Analyses performed with this antibody should be paralleled by positive and negative controls. If unexpected results are obtained which cannot be attributed to differences in laboratory procedures, please contact us. This product may contain hazardous ingredients. Please refer to the Safety Data Sheets (SDS) for additional information and proper handling procedures. Dispose product remainders

according to local regulations. This datasheet is as accurate as reasonably achievable, but Exalpha Biologicals accepts no liability for any inaccuracies or omissions in this information.

## References

1. Ramaekers, F., Huysmans, A., Schaart, G., Moesker, O., and Vooijs, P. (1987). Tissue distribution of keratin 7 as monitored by a monoclonal antibody, *Exp Cell Res* 170, 235-49.
2. Smedts, F., Ramaekers, F., Robben, H., Pruszczynski, M., van Muijen, G., Lane, B., Leigh, I., and Vooijs, P. (1990). Changing patterns of keratin expression during progression of cervical intraepithelial neoplasia, *Am J Pathol* 136, 657-68.
3. Schaafsma, H. E., Ramaekers, F. C., van Muijen, G. N., Lane, E. B., Leigh, I. M., Robben, H., Huijsmans, A., Ooms, E. C., and Ruiter, D. J. (1990). Distribution of cytokeratin polypeptides in human transitional cell carcinomas, with special emphasis on changing expression patterns during tumor progression, *Am J Pathol* 136, 329-43.
4. Ramaekers, F., van Niekerk, C., Poels, L., Schaafsma, E., Huijsmans, A., Robben, H., Schaart, G., and Vooijs, P. (1990). Use of monoclonal antibodies to keratin 7 in the differential diagnosis of adenocarcinomas, *Am J Pathol* 136, 641-55.
5. Bauwens, L. J., De Groot, J. C., Ramaekers, F. C., Veldman, J. E., and Huizing, E. H. (1992). Expression of intermediate filament proteins in the adult Human vestibular labyrinth, *Ann Otol Rhinol Laryngol* 101, 479-86.
6. Smedts, F., Ramaekers, F., Troyanovsky, S., Pruszczynski, M., Link, M., Lane, B., Leigh, I., Schijf, C., and Vooijs, P. (1992). Keratin expression in cervical cancer, *Am J Pathol* 141, 497-511.
7. Ivanyi, D., Minke, J. M., Hageman, C., Groeneveld, E., and van Doornewaard, G. (1992). Patterns of expression of feline cytokeratins in healthy epithelia and mammary carcinoma cells, *Am J Vet Res* 53, 304-14.
8. Glass C, Fuchs E. Isolation, sequence, and differential expression of a human K7 gene in simple epithelial cells. *J Cell Biol.* 1988 Oct;107(4):1337-50. doi: 10.1083/jcb.107.4.1337.
9. Zhou Q, Toivola DM, Feng N, Greenberg HB, Franke WW, Omary MB. Keratin 20 helps maintain intermediate filament organization in intestinal epithelia. *Mol Biol Cell.* 2003 Jul;14(7):2959-71. doi: 10.1091/mbc.e03-02-0059

## CE Mark

CE

## Protein Reference(s)

*Database Name:* UniProt

*Accession Number:* P08729

## Safety Datasheet(s) for this product:

NM\_Sodium Azide



Figure 1 - Immunohistochemistry of MUB0315P (RCK105) on frozen sections of swine liver showing positive staining in the epithelial cells lining the bile ducts and no reactivity in hepatocytes or connective tissue. Dilution 1:500

---



Figure 2 - Immunohistochemistry of MUB0315P (RCK105) on frozen sections of swine liver showing positive staining in the epithelial cells lining the bile ducts and no reactivity in hepatocytes or connective tissue. Dilution 1:500

---



Figure 3 - Immunohistochemistry of MUB0315P (RCK105) on frozen sections of swine liver showing positive staining in the epithelial cells lining the bile ducts and no reactivity in hepatocytes or connective tissue. Dilution 1:500

---



Figure 4. - Immunohistochemistry of MUB0315P (RCK105) on frozen sections of human liver showing positive staining in the epithelial cells lining the bile ducts and no reactivity in hepatocytes or connective tissue. Dilution 1:500

---



Figure 5 - Immunohistochemistry of MUB0315P (RCK105) on frozen section of human urinary bladder showing positive staining in the epithelial cells and no reactivity in connective tissue. Dilution 1:500

---



Figure 6. Immunohistochemistry of MUB0315P (RCK105) on frozen sections of swine liver showing positive staining in the epithelial cells lining the bile ducts and no reactivity in hepatocytes or connective tissue. Dilution 1:500

---



Figure 7 - Immunohistochemistry of MUB0315P (RCK105) on frozen sections of swine liver showing positive staining in the epithelial cells lining the bile ducts and no reactivity in hepatocytes or connective tissue. Dilution 1:500

---



Figure 8 - Immunohistochemistry of MUB0315P (RCK105) on formalin fixed, paraffin embedded sections of human liver showing positive staining in the epithelial cells lining the bile ducts and no reactivity in hepatocytes or connective tissue. Dilution 1:100