



**SZABO
SCANDIC**

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic



Emi1 (m): 293T Lysate: sc-120019

BACKGROUND

Emi1 (for early mitotic inhibitor) regulates mitosis by inhibiting the anaphase promoting complex/cyclosome (APC). Emi1 is a conserved F box protein containing a zinc binding region essential for APC inhibition. The Emi1 protein functions to promote cyclin A accumulation and S phase entry in somatic cells by inhibiting the APC complex. At the G₁-S transition, Emi1 is transcriptionally induced by the E2F transcription factor. Emi1 overexpression accelerates S phase entry and can override a G₁ block caused by overexpression of Cdh1 or the E2F-inhibitor p105 retinoblastoma protein (pRb). Depleting cells of Emi1 through RNA interference prevents accumulation of cyclin A and inhibits S phase entry. Emi1 is required to arrest unfertilized eggs at metaphase of meiosis II and may be the long-sought mediator of CSF activity. Human Emi1 is similar to *Xenopus laevis* Emi1, which inhibits the APC (Cdc20) ubiquitination complex to allow accumulation of cyclin B.

REFERENCES

1. Cenciarelli, C., Chiaur, D.S., Guardavaccaro, D., Parks, W., Vidal, M. and Pagano, M. 1999. Identification of a family of human F-box proteins. *Curr. Biol.* 9: 1177-1179.
2. Chiaur, D.S., Murthy, S., Cenciarelli, C., Parks, W., Loda, M., Inghirami, G., Demetrick, D. and Pagano, M. 2000. Five human genes encoding F-box proteins: chromosome mapping and analysis in human tumors. *Cytogenet. Cell Genet.* 88: 255-258.
3. Reimann, J.D., Freed, E., Hsu, J.Y., Kramer, E.R., Peters, J.M. and Jackson, P.K. 2001. Emi1 is a mitotic regulator that interacts with Cdc20 and inhibits the anaphase promoting complex. *Cell* 105: 645-655.
4. Reimann, J.D. and Jackson, P.K. 2002. Emi1 is required for cytostatic factor arrest in vertebrate eggs. *Nature* 416: 850-854.
5. Jin, J., Cardozo, T., Lovering, R.C., Elledge, S.J., Pagano, M. and Harper, J.W. 2004. Systematic analysis and nomenclature of mammalian F-box proteins. *Genes Dev.* 18: 2573-2580.
6. Eldridge, A.G., Loktev, A.V., Hansen, D.V., Verschuren, E.W., Reimann, J.D. and Jackson, P.K. 2006. The Evi-5 oncogene regulates cyclin accumulation by stabilizing the anaphase-promoting complex inhibitor Emi1. *Cell* 124: 367-380.
7. Lee, H., Lee, D.J., Oh, S.P., Park, H.D., Nam, H.H., Kim, J.M. and Lim, D.S. 2006. Mouse Emi1 has an essential function in mitotic progression during early embryogenesis. *Mol. Cell. Biol.* 26: 5373-5381.
8. Lehman, N.L., Verschuren, E.W., Hsu, J.Y., Cherry, A.M. and Jackson, P.K. 2006. Overexpression of the anaphase promoting complex/cycosome inhibitor Emi1 leads to tetraploidy and genomic instability of p53-deficient cells. *Cell Cycle* 5: 1569-1573.
9. Miller, J.J., Summers, M.K., Hansen, D.V., Nachury, M.V., Lehman, N.L., Loktev, A. and Jackson, P.K. 2006. Emi1 stably binds and inhibits the anaphase-promoting complex/cycosome as a pseudosubstrate inhibitor. *Genes Dev.* 20: 2410-2420.

CHROMOSOMAL LOCATION

Genetic locus: Fbxo5 (mouse) mapping to 10 A1.

PRODUCT

Emi1 (m): 293T Lysate represents a lysate of mouse Emi1 transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

APPLICATIONS

Emi1 (m): 293T Lysate is suitable as a Western Blotting positive control for mouse reactive Emi1 antibodies. Recommended use: 10-20 µl per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.