



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

NDUFV3 (m): 293T Lysate: sc-121986

BACKGROUND

Located in the mitochondrial inner membrane, mitochondrial complex I is the first and largest enzyme in the electron transport chain of oxidative phosphorylation. By oxidizing NADH that is produced in the Krebs cycle, this complex utilizes the two electrons to reduce ubiquinone to ubiquinol, thereby initiating the passage of electrons to successive complexes and ultimately leading to the reduction of oxygen to water. Mitochondrial complex I consists of over 40 subunits and is of considerable clinical interest since defects in any one of the subunits can lead to various myopathies and neuropathies. As a subunit of mitochondrial complex I, NDUFV3 (NADH dehydrogenase [ubiquinone] flavo-protein 3), also designated NADH-ubiquinone oxidoreductase 9 kDa subunit or CI-9kD, is a 108 amino acid protein that is believed to not be involved in catalysis.

REFERENCES

- Pilkington, S.J. and Walker, J.E. 1989. Mitochondrial NADH-ubiquinone reductase: complementary DNA sequences of import precursors of the bovine and human 24 kDa subunit. *Biochemistry* 28: 3257-3264.
- Sled, V.D. and Vinogradov, A.D. 1993. Kinetics of the mitochondrial NADH-ubiquinone oxidoreductase interaction with hexammineruthenium(III). *Biochim. Biophys. Acta* 1141: 262-268.
- de Coo, R.F., Buddiger, P., Smeets, H.J. and van Oost, B.A. 1997. Molecular cloning and characterization of the human mitochondrial NADH:oxidoreductase 10 kDa gene (NDUFV3). *Genomics* 45: 434-437.
- Scanlan, M.J., Gordan, J.D., Williamson, B., Stockert, E., Bander, N.H., Jongeneel, V., Gure, A.O., Jäger, D., Jäger, E., Knuth, A., Chen, Y.T. and Old, L.J. 1999. Antigens recognized by autologous antibody in patients with renal-cell carcinoma. *Int. J. Cancer* 83: 456-464.
- Berry, A., Scott, H.S., Kudoh, J., Talior, I., Korostishevsky, M., Wattenhofer, M., Guipponi, M., Barras, C., Rossier, C., Shibuya, K., Wang, J., Kawasaki, K., Asakawa, S., Minoshima, S., Shimizu, N., Antonarakis, S. and Bonne-Tamir, B. 2000. Refined localization of autosomal recessive nonsyndromic deafness DFNB10 locus using 34 novel microsatellite markers, genomic structure, and exclusion of six known genes in the region. *Genomics* 68: 22-29.
- Zakharova, N.V. and Zharova, T.V. 2002. Kinetic mechanism of mitochondrial NADH:ubiquinone oxidoreductase interaction with nucleotide substrates of the transhydrogenase reaction. *Biochemistry* 67: 1395-1404.
- Kervinen, M., Hinttala, R., Helander, H.M., Kurki, S., Uusimaa, J., Finel, M., Majamaa, K. and Hassinen, I.E. 2006. The MELAS mutations 3946 and 3949 perturb the critical structure in a conserved loop of the ND1 subunit of mitochondrial complex I. *Hum. Mol. Genet.* 15: 2543-2552.
- Ben-Shachar, D. and Karry, R. 2008. Neuroanatomical pattern of mitochondrial complex I pathology varies between schizophrenia, bipolar disorder and major depression. *PLoS ONE* 3: e3676.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

CHROMOSOMAL LOCATION

Genetic locus: Ndufv3 (mouse) mapping to 17 B1.

PRODUCT

NDUFV3 (m): 293T Lysate represents a lysate of mouse NDUFV3 transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

APPLICATIONS

NDUFV3 (m): 293T Lysate is suitable as a Western Blotting positive control for mouse reactive NDUFV3 antibodies. Recommended use: 10-20 µl per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.