



**SZABO
SCANDIC**

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic



p21 (m): 293T Lysate: sc-122305

BACKGROUND

It is now well established that cyclins play a positive role in promoting cell cycle transitions via their ability to associate with and activate their cognate cyclin-dependent kinases (Cdks). Cdk2 associates with cyclins A, D and E, and has been implicated in the control of the G₁ to S phase transition in mammals. A novel Cdk-interacting protein, designated p21, Cip1 or WAF1, has been identified in cyclin A, cyclin D1, cyclin E and Cdk2 immunoprecipitates. p21 is a potent, tight-binding inhibitor of Cdks and can inhibit the phosphorylation of Rb by cyclin A-Cdk2, cyclin E-Cdk2, cyclin D1-Cdk4 and cyclin D2-Cdk4 complexes. Expression of p21 is inducible by wildtype, but not mutant, p53. The mouse homolog of p21 is designated CAP20.

REFERENCES

- Sherr, C.J. 1993. Mammalian G₁ cyclins. *Cell* 73: 1059-1065.
- Harper, J.W., Adami, G.R., Wei, N., Keyomarsi, K. and Elledge, S.J. 1993. The p21 Cdk-interacting protein Cip1 is a potent inhibitor of G₁ cyclin-dependent kinases. *Cell* 75: 805-816.
- El-Deiry, W.S., Tokino, T., Velculescu, V.E., Levy, D.B., Parsons, R., Trent, J.M., Lin, D., Mercer, W.E., Kinzler, K.W. and Vogelstein, B. 1993. WAF1, a potential mediator of p53 tumor suppression. *Cell* 75: 817-825.
- Hunter, T. 1993. Braking the cycle. *Cell* 75: 839-841.
- Kato, J., Matsushime, H., Hiebert, S.W., Ewen, M.E. and Sherr, C.J. 1993. Direct binding of cyclin D to the retinoblastoma gene product and pRb phosphorylation by the cyclin D-dependent kinase Cdk4. *Genes Dev.* 7: 331-342.
- Xiong, Y., Hannon, G.J., Zhang, H., Casso, D., Kobayashi, R. and Beach, D. 1993. p21 is a universal inhibitor of cyclin kinases. *Nature* 366: 701-704.
- Gu, Y., Turck, C.W. and Morgan, D.O. 1993. Inhibition of Cdk2 activity *in vivo* by an associated 20 kDa regulatory subunit. *Nature* 366: 707-710.
- El-Deiry, W.S., Harper, J.W., O'Connor, P.M., Velculescu, V.E., Canman, C.E., Jackman, J., Pietenpol, J.A., Burrell, M., Hill, D.E., Wang, Y., et al. 1994. WAF1/Cip1 is induced in p53-mediated G₁ arrest and apoptosis. *Cancer Res.* 54: 1169-1174.

CHROMOSOMAL LOCATION

Genetic locus: Cdkn1a (mouse) mapping to 17 A3.3.

PRODUCT

p21 (m): 293T Lysate represents a lysate of mouse p21 transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

APPLICATIONS

p21 (m): 293T Lysate is suitable as a Western Blotting positive control for mouse reactive p21 antibodies. Recommended use: 10-20 µl per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.


p21 (B-2): sc-271532 is recommended as a positive control antibody for Western Blot analysis of enhanced mouse p21 expression in p21 transfected 293T cells (starting dilution 1:100, dilution range 1:100-1:1,000).

RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended:

- Western Blotting: use m-IgG_κ BP-HRP: sc-516102 or m-IgG_κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048.

DATA



p21 (B-2): sc-271532. Western blot analysis of p21 expression in non-transfected: sc-117752 (**A**) and mouse p21 transfected: sc-122305 (**B**) 293T whole cell lysates.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.