



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# CYP26C1 siRNA (m): sc-142666

## BACKGROUND

Cytochrome P450 proteins are heme-thiolate monooxygenases that mediate NADPH-dependent electron transport and function to oxidize a variety of structurally unrelated compounds, including steroids, fatty acids and xenobiotics. Specifically, cytochrome P450s are responsible for metabolizing arachidonic acid to hydroxyeicosatetraenoic acid (a regulator of blood pressure) and epoxyeicosatrienoic acid (a molecule involved in signaling events). CYP26C1 (cytochrome P450, family 26, subfamily C, polypeptide 1) is a 522 amino acid single-pass membrane protein expressed in most tissues at very low levels. Belonging to the cytochrome P450 family, CYP26C1 uses heme groups as cofactors and is induced by retinoic acid. CYP26C1 participates in the retinoic metabolism and acts on retinoids, including all-*trans*-retinoic acid (RA) and its stereoisomer, 9-*cis*-RA.

## REFERENCES

1. Ozpolat, B., Mehta, K., Tari, A.M. and Lopez-Berestein, G. 2002. all-*trans*-Retinoic acid-induced expression and regulation of retinoic acid 4-hydroxylase (CYP26) in human promyelocytic leukemia. *Am. J. Hematol.* 70: 39-47.
2. Ozpolat, B. and Lopez-Berestein, G. 2002. Liposomal-all-*trans*-retinoic acid in treatment of acute promyelocytic leukemia. *Leuk. Lymphoma* 43: 933-941.
3. Hercule, H.C., Wang, M.H. and Oyekan, A.O. 2003. Contribution of cytochrome P450 4A isoforms to renal functional response to inhibition of nitric oxide production in the rat. *J. Physiol.* 551: 971-979.
4. Taimi, M., Helvig, C., Wisniewski, J., Ramshaw, H., White, J., Amad, M., Korczak, B. and Petkovich, M. 2004. A novel human cytochrome P450, CYP26C1, involved in metabolism of 9-*cis* and all-*trans* isomers of retinoic acid. *J. Biol. Chem.* 279: 77-85.
5. Nelson, D.R., Zeldin, D.C., Hoffman, S.M., Maltais, L.J., Wain, H.M. and Nebert, D.W. 2004. Comparison of cytochrome P450 (CYP) genes from the mouse and human genomes, including nomenclature recommendations for genes, pseudogenes and alternative-splice variants. *Pharmacogenetics* 14: 1-18.
6. Idres, N., Marill, J. and Chabot, G.G. 2005. Regulation of CYP26A1 expression by selective RAR and RXR agonists in human NB4 promyelocytic leukemia cells. *Biochem. Pharmacol.* 69: 1595-1601.
7. Zhang, F., Chen, C.L., Qian, J.Q., Yan, J.T., Cianflone, K., Xiao, X. and Wang, D.W. 2005. Long-term modifications of blood pressure in normotensive and spontaneously hypertensive rats by gene delivery of rAAV-mediated cytochrome P450 arachidonic acid hydroxylase. *Cell Res.* 15: 717-724.
8. Ozpolat, B., Mehta, K. and Lopez-Berestein, G. 2005. Regulation of a highly specific retinoic acid-4-hydroxylase (CYP26A1) enzyme and all-*trans*-retinoic acid metabolism in human intestinal, liver, endothelial, and acute promyelocytic leukemia cells. *Leuk. Lymphoma* 46: 1497-1506.
9. Kuczeriszka, M., Badzyńska, B. and Kompanowska-Jezińska, E. 2006. Cytochrome P-450 monooxygenases in control of renal haemodynamics and arterial pressure in anaesthetized rats. *J. Physiol. Pharmacol.* 57: 179-185.

## CHROMOSOMAL LOCATION

Genetic locus: Cyp26c1 (mouse) mapping to 19 C2.

## PRODUCT

CYP26C1 siRNA (m) is a pool of 2 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see CYP26C1 shRNA Plasmid (m): sc-142666-SH and CYP26C1 shRNA (m) Lentiviral Particles: sc-142666-V as alternate gene silencing products.

For independent verification of CYP26C1 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-142666A and sc-142666B.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

CYP26C1 siRNA (m) is recommended for the inhibition of CYP26C1 expression in mouse cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor CYP26C1 gene expression knockdown using RT-PCR Primer: CYP26C1 (m)-PR: sc-142666-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.