

# Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

## Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

#### SANTA CRUZ BIOTECHNOLOGY, INC.

## PGRP-I $\alpha$ siRNA (m): sc-152195



#### BACKGROUND

Peptidoglycan recognition proteins (PGRPs) are molecules that recognize peptidoglycan, a large component in bacterial cell walls. In insects, PGRPs activate antimicrobial pathways, and in mammals PGRPs function as antibacterial neutrophil proteins. PGRP-L halts bacterial growth by acting as an alanine amidase, an enzyme that hydrolyzes the amide bond of bacterial peptidoglycan. PGRP-I $\alpha$  and PGRP-I $\beta$  are also members of the PGRP family that help to recognize bacteria by binding to peptidoglycan and Gram-positive bacteria, but they do not have amidase activity. These two PGRPs are expressed in the esophagus and, to a lesser extent, in the tonsils and thymus. PGRP-I $\alpha$  and PGRP-I $\beta$  are transmembrane proteins of 341 and 373 amino acids, respectively, and they have have at least three highly conserved C-terminal PGRP domains either in the extracellular or in the cytoplasmic (or in both) regions.

#### REFERENCES

- Liu, C., Xu, Z., Gupta, D. and Dziarski, R. 2001. Peptidoglycan recognition proteins: a novel family of four human innate immunity pattern recognition molecules. J. Biol. Chem. 276: 34686-34694.
- Wang, Z.M., Li, X., Cocklin, R.R., Wang, M., Wang, M., Fukase, K., Inamura, S., Kusumoto, S., Gupta, D. and Dziarski, R. 2003. Human peptidoglycan recognition protein-L is an N-acetylmuramoyl-L-alanine amidase. J. Biol. Chem. 278: 49044-49052.
- Guan, R., Malchiodi, E.L., Wang, Q., Schuck, P. and Mariuzza, R.A. 2004. Crystal structure of the C-terminal peptidoglycan-binding domain of human peptidoglycan recognition protein Iα. J. Biol. Chem. 279: 31873-31882.
- Natori, S. 2004. Overview: Innate immunity and peptideglycan recognition protein. Tanpakushitsu Kakusan Koso 49: 1156-1160.
- 5. Fournier, B. and Philpott, D.J. 2005. Recognition of *Staphylococcus aureus* by the innate immune system. Clin. Microbiol. Rev. 18: 521-540.
- 6. Kumar, S., Roychowdhury, A., Ember, B., Wang, Q., Guan, R., Mariuzza, R.A. and Boons, G.J. 2005. Selective recognition of synthetic lysine and mesodiaminopimelic acid-type peptidoglycan fragments by human peptidoglycan recognition proteins I{ $\alpha$ } and S. J. Biol. Chem. 280: 37005-37012.
- Uehara, A., Sugawara, Y., Kurata, S., Fujimoto, Y., Fukase, K., Kusumoto, S., Satta, Y., Sasano, T., Sugawara, S. and Takada, H. 2005. Chemically synthesized pathogen-associated molecular patterns increase the expression of peptidoglycan recognition proteins via toll-like receptors, NOD1 and NOD2 in human oral epithelial cells. Cell. Microbiol. 7: 675-686.
- Wang, H., Gupta, D., Li, X. and Dziarski, R. 2005. Peptidoglycan recognition protein 2 (N-acetylmuramoyl-L-Ala amidase) is induced in keratinocytes by bacteria through the p38 kinase pathway. Infect. Immun. 73: 7216-7225.
- 9. Guan, R., Brown, P.H., Swaminathan, C.P., Roychowdhury, A., Boons, G.J. and Mariuzza, R.A. 2006. Crystal structure of human peptidoglycan recognition protein I  $\alpha$  bound to a muramyl pentapeptide from gram-positive bacteria. Protein Sci. 15: 1199-1206.

#### CHROMOSOMAL LOCATION

Genetic locus: Pglyrp3 (mouse) mapping to 3 F1.

#### PRODUCT

PGRP-I $\alpha$  siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see PGRP-I $\alpha$  shRNA Plasmid (m): sc-152195-SH and PGRP-I $\alpha$  shRNA (m) Lentiviral Particles: sc-152195-V as alternate gene silencing products.

For independent verification of PGRP-I $\alpha$  (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-152195A, sc-152195B and sc-152195C.

#### STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNAse-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNAse-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

#### **APPLICATIONS**

PGRP-I  $\alpha$  siRNA (m) is recommended for the inhibition of PGRP-I  $\alpha$  expression in mouse cells.

#### SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

#### **RT-PCR REAGENTS**

Semi-quantitative RT-PCR may be performed to monitor PGRP-I $\alpha$  gene expression knockdown using RT-PCR Primer: PGRP-I $\alpha$  (m)-PR: sc-152195-PR (20  $\mu$ I). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

#### **RESEARCH USE**

For research use only, not for use in diagnostic procedures.

#### PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.