



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Aluminium cryo tubes + caps, 15 ml, 29mmx 30mm

Reagents

Instructions for use

Product name	Aluminium cryo tubes + caps, 15 ml, 29mmx 30mm
Intended Use	Storage of tissue samples. The aluminium cryo vials with screw caps can be used for storage of tissue samples. The cryo tubes are suitable for storage of tissue samples in liquid nitrogen. Suitable for deepfreezing till -196°C. Autoclavable at 121°C . Available in 3ml (PA6003) and 15ml (PA6015) (sets of 100 pieces).
Applications	IHC-P, IHC-Fr
Summary and explanation	-

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Aluminium cryo tubes + caps, 15 ml, 29mmx 30mm

Reagents

Instructions for use

Principle of method

-

Reagents provided

100 tubes with caps

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com

FOR-029 22-12-2021

Revision date

02-08-2022

Product datasheet

PA6015

MONOSAN[®]

Aluminium cryo tubes + caps, 15 ml, 29mmx 30mm

Reagents

Instructions for use

Storage and handling Storage of cryo tubes in liquid nitrogen: 1) Place the tissue in the cryo tube, 2) Close the cryotube tightly and place the cryotube in the liquid nitrogen for 1 minute, 3) The cryo vials can be stored at -80°C. Use of cryo tube in combination with isopentane: 1) Place a cup with isopentane in the liquid nitrogen; after approximately 2 minutes the isopentane clots and has a temperature of -160°C. 2) Place the tissue sample or biopsy in the isopentane for 30 seconds, 3) Place the tissue in the cryo vials and store at -80°C.

Reagent preparation Make sure that the cap and the vial have the same temperature because the material could slightly shrink or expand under influence of excessive temperature differences.

Procedure -

Expected results Long term storage of tissue

Trouble shooting

-

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com

FOR-029 22-12-2021

Revision date
02-08-2022

Product datasheet

PA6015

MONOSAN[®]

Aluminium cryo tubes + caps, 15 ml, 29mmx 30mm

Reagents

Instructions for use

Quality control -

Performance -

Limitations of procedure -

Precautions

Wear a face mask at all times, Check if the screw treads and the lid are okay, Screw the lid on tightly. Attention: AVOID EXCESSIVE TEMPERATURE CHANGES, the cap could explode from the cryo tube.

References

1. -
2. -
3. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES