



**SZABO
SCANDIC**

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic



CD23, B-Cell, (Clone MHM6)

Description:

Species: Mouse, Monoclonal.

Immunogen: EBV-transformed B lymphoblastoid cells.

Clone: MHM6.

Isotype: IgG1, Kappa.

Format: This antibody has been pretitered and quality controlled to work on formalin-fixed paraffin-embedded as well as acetone fixed cryostat tissue sections. No further titration is required.

Specificity: In lymphoid tissue, the antibody labels a variable proportion of mantle zone lymphocytes and a fraction of dendritic reticulum cells. In lymph nodes with progressive germinal centre transformation, the number and staining intensity of small follicular lymphocytes are increased. Essentially no staining is seen on splenic marginal zone lymphocytes.

The CD23 antigen is expressed on neoplastic cells from cases of B-Cell chronic lymphocytic leukemia and some cases of centroblastic/centrocytic lymphoma, but not in other types of lymphoid neoplasms.

A large number of normal non-lymphoid human tissues do not react with CD23, (Clone MHM6) including: central nervous system, endocrine glands, gastrointestinal, respiratory, urinary, reproductive system, skin, bone, cartilage, connective tissue, muscle, peripheral nerves, vessels, mesothelium, synovium, placenta and umbilical cord.

Uses/Limitations: Immunohistochemistry
For Research Use Only.
Do not use past expiration date.

Storage: 2-8° Centigrade.

Procedure: We suggest an incubation period of 30-60 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the user. For formalin-fixed paraffin-embedded tissue pretreatment with Citrate Plus (catalog # CPL500) is **required**.

Precautions: Contains Sodium Azide as a preservative.
Do not pipette by mouth.

References:

- 1) Kukutani et al. White Cell Differentiation Antigens, 419-22, 1987.
- 2) Sarfati et al. White Cell Differentiation Antigens, 530-3, 1995.
- 3) Thorely-Lawson et al. J Immunol 134: 3007-12, 1985.
- 4) Defrance et al. J Exp Med, 179: 135-43, 1994
- 5) Defrance et al. J Exp Med 179: 135-43, 1994

Instructions For Use A00081-IFU

Rev. Date: 10/15/03 **Revision: 3** Page 2 of 2

P.O. Box 3286 - Logan, Utah 84323, U.S.A. - Tel. (800) 729-8350 - Fax (435) 755-0015 - www.scytek.com

- 6) Aubry et al. White Cell Differentiation Antigens, 417-19, 1987.
- 7) Hardie et al. Eur J Immunol 23: 997-1004, 1993.
- 8) Kikutani et al. J Exp Med 164: 1455-69, 1986
- 9) Rowe et al. Int J Cancer 29: 373-81, 1982
- 10) Pallesen et al. White Cell Differentiation Antigens, 383-6, 1987.