



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## DAB Enhancer Solution

**Description:** 3,3'Diaminobenzidine (DAB) is a widely used chromogen for immunohistochemical staining and immunoblotting. When in the presence of peroxidase enzyme, DAB produces a brown precipitate that is insoluble in alcohol. This solution has been developed to enhance DAB staining. The DAB will stain darker red/brown.

**Availability/Contents:**

<u>Item #</u>	<u>Volume</u>
ACM015	15ml Dropper
ACM125	125ml
ACM500	500ml
ACM999	1000ml

**Storage:** Room Temperature.

**Procedure:**

1. Prior to DAB incubation, add 20  $\mu$ l of DAB Enhancer to 1 ml of DAB ready-to-use reagent.
2. Stain with DAB reagent as usual.
3. Rinse tissue, counterstain, and coverslip.