



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Phosphate Buffered Saline Plus Tween (10x Solution)

Description: Phosphate Buffered Saline Plus Tween has been formulated in a 10X concentration for use with immunolabeling techniques as a wash buffer. The need to weigh and mix individual components is eliminated and replaced with a simple dilution.

Contents: Phosphate buffered saline in reagent grade water with diluted concentrations of Sodium Chloride, NaCl: 0.15M; Sodium Phosphate, Na₂HPO₄: 8mM; Potassium Chloride, KCl: 3mM; Potassium Phosphate, Monobasic, KH₂PO₄: 2mM, Tween, 0.05 %. Final pH of concentrated buffer is 7.4±0.05, diluted buffer is 7.9±0.2. This product is filtered to 0.2 microns.

Availability:

<u>Item #</u>	<u>Volume</u>	<u>DILUTED VOLUME</u>
PBE999	1 liter	10 liters
PBE010	10 liters	100 liters

Storage: Store Phosphate Buffered Saline at room temperature. Product is stable for 24 months from date of manufacture.

Diluting Procedure:

1. Pour Phosphate Buffered Saline plus Tween (10x) in mixing flask and add water to final volume. Stir briefly.

Precautions: Do not pipet reagent by mouth.