



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Horseradish Peroxidase Avidin D (Av-HRP), Concentrate, for IHC

A-2004-5

[Product Images](#)



Short Description

Horseradish Peroxidase Avidin D (HRP Avidin D) is produced by our own coupling procedure that preserves the high specific activity of the peroxidase. Using biotinylated primary or secondary intermediates and peroxidase-labeled Avidin D, antigens can be localized in histological sections, cytospin preparations, or smears.

Features:

- Avidin and streptavidin reagents are powerful tools to detect or purify biotinylated proteins, nucleic acids, and other macromolecules.
- Horseradish Peroxidase Avidin D is offered as a concentrate (5 mg/ml) or in a ready-to-use (R.T.U.) stabilized solution (100ml, 1 µg/ml) in a bottle fitted with a drop dispenser top.

Additional Information

Unit Size	5 mg
Applications	Immunohistochemistry / Immunocytochemistry, In situ hybridization
Recommended Usage	For diluting this product, we recommend a buffered saline solution, approximately pH 7.5. The recommended concentration range for use is 2-5 µg/ml.
Solution	0.05 M NaHCO ₃ , pH 8.2, preservative
Recommended Storage	2-8 °C DO NOT FREEZE
Concentration	5 mg active conjugate/ml
Detection Enzyme(s)	Peroxidase
Conjugate	Peroxidase
Format	Concentrate

