



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



---

# Streptavidin, Agarose

## SA-5010-2

Product Images

---



## Short Description

---

Agarose Streptavidin (SA-5010) can be used to separate biotinylated macromolecules from unbiotinylated materials or for solid phase binding assays.

### Features:

- Biotin binding activity of the streptavidin optimally preserved through proprietary coupling procedure
- Low non-specific binding activity of gel
- Protein concentration is 1 mg streptavidin per ml of settled agarose beads
- Supplied as a 1:1 suspension in buffer
- Binding capacity is approximately 1 mole biotinylated HRP per mole streptavidin

## Additional Information

---

Unit Size	2 ml
Applications	Affinity Chromatography
Concentration	1 mg/ml of settled gel
Recommended Storage	2-8 °C, DO NOT FREEZE
Solution	10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide.
Recommended Usage	Wash gel thoroughly with buffer before use.
Matrix Conjugate	Streptavidin

