



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Levamisole Solution

Cat. No. SP-5000

Storage Refrigerate 2-8 °C

Description Levamisole will inhibit most forms of alkaline phosphatase other than the intestinal isoenzyme. Since the intestinal form of the enzyme is used as the marker enzyme in VECTASTAIN® ABC-AP Kits, levamisole can be used to inhibit endogenous alkaline phosphatase in most cases. Because the inhibition produced by levamisole is reversible, the inhibitor must be added to the substrate solution.

Components

Product Name	Volume
Levamisole Solution	18 ml

Guidelines for Use

When used in conjunction with Vector Laboratories' alkaline phosphatase substrate kits, add one drop of Levamisole Solution to 5 ml of alkaline phosphatase substrate buffer and mix well. Then add substrate reagents and continue as outlined in the substrate kit instructions.

Detailed product listings, specifications, protocols and additional information is available on our website: vectorlabs.com