

official distributor

SZABO-SCANDIC Handels GmbH
Quellenstraße 110, A-1100 Wien
T. +43(0)1 489 3961-0
F. +43(0)1 489 3961-7
mail@szabo-scandic.com
www.szabo-scandic.com



SZABO
SCANDIC



*Just add DNA and Go
for MultiPlex PCR*





Founded in 1987, Clonit srl was one of the first European companies operating in the field of biotechnology. Since from the beginning the company has worked to develop innovative and reliable analytical methods in the field of In Vitro diagnostics using molecular biology techniques, with particular attention to the detection of viral, bacteria and protozoa infections, genetic mutations and oncological diseases.

The ready to use reagents are suitable for quantitative or qualitative application on most common real time analyzers, as the Lifetechnologies 7000 serie, quantstudio serie and StepOne, the CFX 96 from Biorad, the Rotor-GeneQ from Qiagen, the Roche LightCycler 480, the Stratagene from Agilent technologies, the Versant kPCR from Siemens and MIC from BMS. The kits include all the reagents needed to perform the assays, including the standards to built the calibration curve and positive and negative controls.

The standardization of internal quality control system and cross-checks in outer reference labs have made Clonit able to assure the highest standards of care and quality (ISO 9001 and ISO 13485 certified).

Clonit offers to its customers an unrivalled experience in molecular biology.





Products

Clonit'ngo is the Ready to Use solution designed for the diagnosis of infectious diseases caused by different pathogens in human samples.

Clonit'ngo combines the power of Multiplex PCR Detection, with a Ready to Use and Room Temperature Storage Format, allowing cost control by single test solution approach.

Key Features

- › **Same thermal protocol for all Clonit'ngo Kits**
- › **Lyophilised product. Just Add Buffer, DNA..'ngo**
- › **Ready & Easy-to-use Kits**
- › **Transport and storage at room temperature**
- › **Shelf-life: 24 months**
- › **1 well / 1 Test**
- › **High sensibility, specificity and reproducibility**
- › **Validated according to ISO 13485 and CE marked**



MonoPlex Sexual Health

CAT#	Product Description	TESTS
CLNG-96-56	Clonit'ngo Streptococcus B	96
CLNG-96-61	Clonit'ngo Treponema pallidum	96



MultiPlex Sexual Health

CAT#	Product Description	TESTS
CLNG-48-00	Clonit'ngo Sexually Transmitted Diseases Top 7 Neisseria gonorrhoeae, Chlamydia trachomatis, Mycoplasma genitalium, Trichomonas vaginalis, Ureaplasma urealyticum, Ureaplasma parvum and/or Mycoplasma hominis	48
CLNG-96-57	Clonit'ngo Herpes Virus 1, Herpes Virus 2 & Treponema pallidum	96
CLNG-96-58	Clonit'ngo H. ducreyi + C. trachomatis	96
CLNG-96-59	Clonit'ngo C. albicans, G. vaginalis & T. vaginalis	96
CLNG-96-60	Clonit'ngo Papilloma Virus 16 + 18	96



Clonit'ngo MultiPlex Gastrointestinal

CAT#	Product Description	TESTS
CLNG-96-18	<i>Clonit'ngo</i> Aeromonas + Yersinia	96
CLNG-96-16	<i>Clonit'ngo</i> Salmonella, Campylobacter & Shigella/EIEC	96
CLNG-96-17	<i>Clonit'ngo</i> Salmonella, Campylobacter & Yersinia	96
CLNG-96-19	<i>Clonit'ngo</i> E. coli, EPEC & EIEC	96
CLNG-96-20	<i>Clonit'ngo</i> E. coli, ETEC & EIEC	96
CLNG-48-21	<i>Clonit'ngo</i> E. coli Typing	48
CLNG-96-07	<i>Clonit'ngo</i> Norovirus GI + GII	96
CLNG-96-26	<i>Clonit'ngo</i> Cryptosporidium, Giardia & Entamoeba histolytica	96
CLNG-96-28	<i>Clonit'ngo</i> Blastocystis hominis + Dientamoeba fragilis	96
CLNG-96-64	<i>Clonit'ngo</i> H. pylori + Clarithromycin resistance	96 new



Clonit'ngo MonoPlex Gastrointestinal

CAT#	Product Description	TESTS
CLNG-96-01	Clonit'ngo Rotavirus	96
CLNG-96-02	Clonit'ngo Adenovirus	96
CLNG-96-03	Clonit'ngo Astrovirus	96
CLNG-96-06	Clonit'ngo Sapovirus	96
CLNG-96-08	Clonit'ngo Helicobacter pylori	96
CLNG-96-09	Clonit'ngo Clostridium difficile	96
CLNG-96-10	Clonit'ngo Clostridium difficile TOX B	96
CLNG-96-11	Clonit'ngo Clostridium difficile TOX A/B	96
CLNG-96-12	Clonit'ngo Campylobacter	96
CLNG-96-13	Clonit'ngo Salmonella	96
CLNG-96-14	Clonit'ngo Shigella/EIEC	96
CLNG-96-15	Clonit'ngo Yersinia enterocolitica	96
CLNG-96-22	Clonit'ngo Cryptosporidium	96
CLNG-96-23	Clonit'ngo Giardia lamblia	96
CLNG-96-24	Clonit'ngo Entamoeba histolytica	96
CLNG-96-25	Clonit'ngo Entamoeba dispar	96
CLNG-96-27	Clonit'ngo Dientamoeba fragilis	96



Clonit'ngo Vector MonoPlex

CAT#	Product Description	TESTS
CLNG-96-70	Clonit'ngo Dengue Virus	96
CLNG-96-75	Clonit'ngo Mayaro Virus	96
CLNG-96-76	Clonit'ngo Yellow Fever	96
CLNG-96-79	Clonit'ngo Trypanosoma cruzi (Chagas)	96 new
CLNG-96-80	Clonit'ngo Crimean-Congo Hemorrhagic Fever Virus	96
CLNG-96-81	Clonit'ngo Japanese Encephalitis Virus	96



Clonit'ngo Vector MultiPlex

CAT#	Product Description	TESTS
CLNG-96-73	Clonit'ngo Zika, Dengue & Chikungunya	96
CLNG-96-78	Clonit'ngo Borrelia, Anaplasma & Coxiella	96
CLNG-48-71	Clonit'ngo Dengue Serotyping kit	48
CLNG-96-77	Clonit'ngo Tick Borne Diseases (7 pathogens) TBEV, Rickettsia spp., Babesia (microti & divergens), Ehrlichia (chafeensis & muris), Borrelia (burgdorferi sensu lato & miyamotoi & hermsii), Anaplasma phagocitophylum e Coxiella burnetii	48 new



MonoPlex Respiratory

CAT#	Product Description	TESTS
CLNG-96-29	<i>Clonit'ngo</i> Influenza A	96
CLNG-96-30	<i>Clonit'ngo</i> Influenza B	96
CLNG-96-32	<i>Clonit'ngo</i> Influenza H1N1	96
CLNG-96-40	<i>Clonit'ngo</i> Parainfluenza 1	96
CLNG-96-41	<i>Clonit'ngo</i> Parainfluenza 2	96
CLNG-96-42	<i>Clonit'ngo</i> Parainfluenza 3	96
CLNG-96-43	<i>Clonit'ngo</i> Parainfluenza 4	96
CLNG-48-44	<i>Clonit'ngo</i> Parainfluenza Typing	48
CLNG-96-45	<i>Clonit'ngo</i> Bocavirus	96
CLNG-96-51	<i>Clonit'ngo</i> Legionella pneumophila	96
CLNG-96-46	<i>Clonit'ngo</i> Human Metapneumovirus	96
CLNG-96-55	<i>Clonit'ngo</i> Pneumocystis jiroveci	96



Clonit'ngo MultiPlex Respiratory

CAT#	Product Description	TESTS
CLNG-96-35	Clonit'ngo RSV A+ RSV B	96
CLNG-96-31	Clonit'ngo Flu A+ Flu B	96
CLNG-96-36	Clonit'ngo Flu A + Flu B +RSV	96
CLNG-96-37	Clonit'ngo Flu Typing I (H1N1 + H3N2)	96
CLNG-96-38	Clonit'ngo Flu Typing II (H1N1, H5N1, H3N2 & H7N9)	96
CLNG-96-52	Clonit'ngo C.pneumoniae, M.pneumoniae & L.pneumophila	96
CLNG-96-47	Clonit'ngo Adenovirus, Metapneumovirus & Bocavirus	96
CLNG-96-54	Clonit'ngo H. Influenzae + S. pneumoniae + M. catarrhalis	96
CLNG-96-48	Clonit'ngo Rhinovirus + Enterovirus	96
CLNG-96-50	Clonit'ngo Coronavirus (229E + NL63 +OC43 & HKU 1)	96
CLNG-48-39	Clonit'ngo Respiratory Viral Panel I (FluA, FluB&RSV)+(H1N1, H5N1, H3N2, H7N9)	48
CLNG-48-49	Clonit'ngo MERS Coronavirus	48
CLNG-96-53	Clonit'ngo Bordetella (B. pertussis, B. parapertussis, B. holmesii)	96




Resistance MultiPlex

CAT#	Product Description	TESTS
CLNG-96-65	Clonit'ngo Vancomycin resistance	96
CLNG-96-64	Clonit'ngo H.pylori + Clarithromycin resistance	96



Meningitis & Transplantation MultiPlex

CAT#	Product Description	TESTS
CLNG-96-62	Clonit'ngo H. Influenzae, N.meningitis & S. pneumoniae	96
CLNG-96-63	Clonit'ngo S. agalactiae, L. monocytogenes & E. coli	96
CLNG-96-66	Clonit'ngo HSV 1, HSV 2 & VZV	96
CLNG-96-67	Clonit'ngo HHV 6,7,8	96
CLNG-96-68	Clonit'ngo BK+ JCV	96



official distributor

SZABO-SCANDIC HandelsgmbH
Quellenstraße 110, A-1100 Wien
T. +43(0)1 489 3961-0
F. +43(0)1 489 3961-7
mail@szabo-scandic.com
www.szabo-scandic.com



SZABO SCANDIC



CLONIT s.r.l.
20139 Milano | Via Quaranta 57 | Italy
tel: +39.02.56814413
fax: +39.02.56814515

info@clonit.it
www.clonit.it