



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7


[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# AlteTube™

## Screw cap tubes in SBS rack format



**0.5ml Screw cap**      **0.75ml Screw cap**      **1.1ml Screw cap**      **1.4ml Screw cap**

	0.5ml Screw cap	0.75ml Screw cap	1.1ml Screw cap	1.4ml Screw cap
Thread type	Internal	Internal	Internal	Internal
Max working volume - ml	0.3	0.49	0.75	0.92
Min storage temperature - °C	-196	-196	-196	-196
Tube height - mm	20	27	36.5	44
Tube height with cap - mm	27	34	43.5	51
Inner diameter - mm	6.8	6.8	6.8	6.8
Outer diameter with cap - mm	8.4	8.4	8.4	8.4
Pitch - mm	9.0	9.0	9.0	9.0
Identification	2D DataMatrix barcode on base			
2D DataMatrix barcode	White on black, 3mm x 3mm, 2 alpha prefix (AL) and 8 numeric Custom prefix available / Black on white available			
Capping options	Screw caps			
Sterility options	AlteClean™ / Sterile (E-beam) / Ethylene oxide (EtO) on request			
Rack format	8 x 12 - SBS (96 tubes)			
Rack material	Polypropylene base with polycarbonate lid			
Rack height with lid - mm	31	37	50	55.5
Packaging options	Bulk packed tubes, no cap – 960/case Bulk packed tubes, precapped – 960/case Racked tubes, no cap – 96/rack, 10 racks/case Racked tubes, precapped – 96/rack, 10 racks/case Bulk screw caps – 960/case Empty racks – 10 racks/case			

Maximum working volumes are for guidance only.

Please request samples for evaluation.

Warning: Caps should not be immersed in Liquid Nitrogen.

AlteTubes should be stored in gas phase Liquid Nitrogen (GPLN) only.

# High quality 2D barcoded tubes to assist sample identification



## SBS format racks

Store 96 tubes in 8 x12 format.

## Rack identification

One unique code in three formats. Embedded 2D barcode on the base of the rack with matching 1D barcode and human readable code.

## AlteClean™

All tubes are manufactured from medical grade virgin polypropylene and certified DNase, RNase and Endotoxin (Pyrogen) free.

## Internal thread screw caps

Manufactured by dual injection of ThermoPlastic Elastomer (TPE) and medical grade polypropylene.

## Exceptional performance

The single piece cap provides a superior fixed seal, eliminating the risk of leakage or cross-threading.

## Capping flexibility

Internal thread AlteTubes are compatible with screw caps and TPE septum caps.



## Permanent 2D barcode

The 2D DataMatrix barcode is laser etched on the black base of the single piece AlteTube.

## 100% quality checks

Every barcode is checked to guarantee uniqueness across the full AlteTube range.



See the full range of 2D barcode readers in the AltemisLab™ product catalogue.

[www.altemislab.com](http://www.altemislab.com)