

Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

PRODUCT INFORMATION



KN-62

Item No. 13318

CAS Registry No.: Formal Name:	127191-97-3 5-isoquinolinesulfonic acid, 4-[(2S)-2-[(5-isoquinolinylsulfonyl) methylamino]-3-oxo-3-(4-phenyl- 1-piperazinyl)propyl]phenyl ester	
MF:	$C_{38}H_{35}N_5O_6S_2$	
FW:	721.9	
Purity:	≥98%	0=s=0
Supplied as:	A crystalline solid	
Storage:	-20°C	
Stability:	≥4 years	

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Laboratory Procedures

KN-62 is supplied as a crystalline solid. A stock solution may be made by dissolving the KN-62 in the solvent of choice, which should be purged with an inert gas. KN-62 is soluble in organic solvents such as ethanol, DMSO, and dimethyl formamide (DMF). The solubility of KN-62 in these solvents is approximately 1 mg/ml in ethanol and 30 mg/ml in DMSO and DMF.

KN-62 is sparingly soluble in aqueous buffers. For maximum solubility in aqueous buffers, KN-62 should first be dissolved in DMF and then diluted with the aqueous buffer of choice. KN-62 has a solubility of approximately 0.3 mg/ml in a 1:2 solution of DMF:PBS (pH 7.2) using this method. We do not recommend storing the aqueous solution for more than one day.

Description

KN-62 is a selective, cell permeable inhibitor of Ca²⁺/calmodulin-dependent kinase type II (CaMKII; IC₅₀ = 900 nM).¹ It does not affect the activity of several other kinases when tested at 10 μ M.² KN-62 also acts as a non-competitive antagonist of the purinergic receptor P2RX7 (IC_{50} = 15 nM).³

References

- 1. Tokunitsu, H., Chijiwa, T., Hagiwara, M., et al. KN-62, 1-[N,O-bis(5-isoquinolinesulfonyl)-N-methyl-L-tyrosyl]-4-phenylpiperazine, a specific Inhibitor of Ca²⁺/calmodulin-dependent protein kinase II. J. Biol. Chem. 265(8), 4315-4320 (1990).
- 2. Davies, S.P., Reddy, H., Caivano, M., et al. Specificity and mechanism of action of some commonly used protein kinase inhibitors. Biochem. J. 351(1), 95-105 (2000).
- 3. Chessell, I.P., Michel, A.D., and Humphrey, P.P.A. Effects of antagonists at the human recombinant P2X₇ receptor. Br. J. Pharmacol. 124(6), 1314-1320 (1998).

WARNING THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFETY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

uyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 11/21/2022

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897 [734] 971-3335 FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.CAYMANCHEM.COM