

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



PRODUCT INFORMATION



Argatroban

Item No. 19639

CAS Registry No.: 74863-84-6 Formal Name: (2R,4R)-1-[(2S)-5-

> [(aminoiminomethyl)amino]-1oxo-2-[[(1,2,3,4-tetrahydro-3methyl-8-quinolinyl)sulfonyl] amino|pentyl]-4-methyl-2piperidinecarboxylic acid

MCI-9038 Synonym: MF: $C_{23}H_{36}N_6O_5S$ FW: 508.6 **Purity:** ≥95%

UV/Vis.: λ_{max} : 212, 259, 331 nm A crystalline solid Supplied as:

-20°C Storage: Stability: ≥2 years

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Laboratory Procedures

Argatroban is supplied as a crystalline solid. A stock solution may be made by dissolving the argatroban in the solvent of choice. Argatroban is soluble in organic solvents such as ethanol, DMSO, and dimethyl formamide, which should be purged with an inert gas. The solubility of argatroban in these solvents is approximately 1, 3, and 20 mg/ml, respectively.

Further dilutions of the stock solution into aqueous buffers or isotonic saline should be made prior to performing biological experiments. Ensure that the residual amount of organic solvent is insignificant, since organic solvents may have physiological effects at low concentrations. Organic solvent-free aqueous solutions of argatroban can be prepared by directly dissolving the crystalline solid in aqueous buffers. The solubility of argatroban in PBS, pH 7.2, is approximately 0.1 mg/ml. We do not recommend storing the aqueous solution for more than one day.

Description

Argatroban is a reversible inhibitor of thrombin (IC $_{50}$ S = 0.01 and 0.09 μ M for the free and clot-bound forms, respectively).1 It is selective for thrombin over factor Xa, trypsin, plasmin, and tissue plasminogen activator (Kis = 0.0045, 53, 0.19, 25.7, and 87.7 µM, respectively). It inhibits thrombin-induced platelet aggregation in washed isolated guinea pig platelets with an IC₅₀ value of 0.077 μM.² Argatroban (3.2 mg/kg, s.c.) decreases thrombus formation in a guinea pig model of iron chloride-induced arterial thrombosis. Formulations containing argatroban have been used in the prevention or treatment of thrombosis in patients with heparin-induced thrombocytopenia.

References

- 1. Park, H.D., Lee, S.H., Kim, T.H., et al. Antithrombotic effects of LB30870, a potent, orally active, selective and direct thrombin inhibitor, and pharmacokinetics of its prodrug. Bioorg. Med. Chem. Lett. 23(17), 4779-4784 (2013).
- 2. Kato, Y., Kita, Y., Hirasawa-Taniyama, Y., et al. Inhibition of arterial thrombosis by a protease-activated receptor 1 antagonist, FR171113, in the guinea pig. Eur. J. Pharmacol. 473(2-3), 163-169 (2003).

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website

Copyright Cayman Chemical Company, 09/08/2020

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897

[734] 971-3335

FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.**CAYMANCHEM**.COM