

Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

PRODUCT INFORMATION



Moxonidine

Item No. 25639

CAS Registry No.:	
Formal Name:	4-chloro-N-(4,5-dihydro-1H-
	imidazol-2-yl)-6-methoxy-2- H H Cl methyl-5-pyrimidinamine \downarrow
Synonym:	LY326869
MF:	$C_0H_{12}CIN_5O$
FW:	241.7 V
Purity:	≥95%
UV/Vis.:	λ _{max} : 222, 256 nm
Supplied as:	A crystalline solid
Storage:	-20°C
Stability:	≥4 years
Information reprocent	the product energifications. Batch energific analytical results are provided on each cortificate of and

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Laboratory Procedures

Moxonidine is supplied as a crystalline solid. A stock solution may be made by dissolving the moxonidine in the solvent of choice, which should be purged with an inert gas. Moxonidine is soluble in organic solvents such as ethanol, DMSO, and dimethyl formamide (DMF). The solubility of moxonidine in ethanol is approximately 10 mg/ml and approximately 20 mg/ml in DMSO and DMF.

Further dilutions of the stock solution into aqueous buffers or isotonic saline should be made prior to performing biological experiments. Ensure that the residual amount of organic solvent is insignificant, since organic solvents may have physiological effects at low concentrations. Organic solvent-free aqueous solutions of moxonidine can be prepared by directly dissolving the crystalline solid in aqueous buffers. The solubility of moxonidine in PBS, pH 7.2, is approximately 2 mg/ml. We do not recommend storing the aqueous solution for more than one day.

Description

Moxonidine is an imidazoline 1 (I_1) receptor agonist ($K_i = 56$ nM in a radioligand binding assay).¹ It is also an agonist of the α_{2A} -adrenergic receptor (α_{2A} -AR; K_i = 150 nM in a radioligand binding assay). It is selective for I₁ and α_{2A} -AR receptors over α_1 -, α_{2B} , and α_{2C} -ARs (K_is = >30,000, 1,000, and 2,000 nM, respectively). Moxonidine decreases mean arterial pressure in spontaneously hypertensive rats when administered at a dose of 0.2 nmol injected into the rostral ventrolateral medulla oblongata (RVLM).² It also reduces blood pressure in cynomolgus monkeys when administered at a dose of 167 μ g/kg.¹

References

- 1. Munk, S.A., Lai, R.K., Burke, J.E., et al. Synthesis and pharmacologic evaluation of 2-endo-amino-3exo-isopropylbicyclo[2.2.1]heptane: A potent imidazoline1 receptor specific agent. J. Med. Chem. 39(6), 1193-1195 (1996).
- 2. Ernsberger, P., Haxhiu, M.A., Graff, L.M., et al. A novel mechanism of action for hypertension control: Moxonidine as a selective I₁-imidazoline agonist. Cardiovasc. Drugs Ther. 8(Suppl 1), 27-41 (1994).

WARNING THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFETY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

uyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 12/07/2022

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897 [734] 971-3335 FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.CAYMANCHEM.COM