

Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

PRODUCT INFORMATION



2-Methylcitric Acid

Item No. 29327

CAS Registry No.:	6061-96-7	
Formal Name:	2-hydroxy-1,2,3-butanetricarboxylic acid	,OH
MF:	$C_7 H_{10} O_7$	
FW:	206.2	
Purity:	≥95%	но он
Supplied as:	A crystalline solid	
Storage:	-20°C	
Stability:	≥2 years	
lu farma atian manufactuta da atiana di atiana Datah ang ifa ang hai atian manufactuta da ang ang ifa da ang a		

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Laboratory Procedures

2-Methylcitric acid is supplied as a crystalline solid. A stock solution may be made by dissolving the 2-methylcitric acid in the solvent of choice, which should be purged with an inert gas. 2-Methylcitric acid is soluble in organic solvents such as ethanol, DMSO, and dimethyl formamide. The solubility of 2-methylcitric acid in these solvents is approximately 10, 30, and 25 mg/ml, respectively.

Further dilutions of the stock solution into aqueous buffers or isotonic saline should be made prior to performing biological experiments. Ensure that the residual amount of organic solvent is insignificant, since organic solvents may have physiological effects at low concentrations. Organic solvent-free aqueous solutions of 2-methylcitric acid can be prepared by directly dissolving the crystalline solid in aqueous buffers. The solubility of 2-methylcitric acid in PBS, pH 7.2, is approximately 10 mg/ml. We do not recommend storing the aqueous solution for more than one day.

Description

2-Methylcitric acid is an endogenous tricarboxylic acid formed by the condensation of propionyl-CoA with oxaloacetic acid by citrate synthase under conditions of propionyl-CoA accumulation.¹ Accumulation of 2-methylcitric acid is associated with cobalamin deficiencies, propionic acidemia, and methylmalonic acidurias.1,2

References

1. Allen, R.H., Stabler, S.P., Savage, D.G., et al. Elevation of 2-methylcitric acid I and II levels in serum, urine, and cerebrospinal fluid of patients with cobalamin deficiency. Metabolism 42(8), 978-988 (1993).

2. Al-Dirbashi, O.Y., McIntosh, N., and Chakraborty, P. Quantification of 2-methylcitric acid in dried blood spots improves newborn screening for propionic and methylmalonic acidemias. J. Med. Screen. 24(2), 58-61 (2017).

WARNING THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFFTY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

uyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 05/12/2020

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897 [734] 971-3335 FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.CAYMANCHEM.COM