

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



PRODUCT INFORMATION



Lidocaine-d₁₀ Item No. 35066

CAS Registry No.: 851528-09-1

Formal Name: 2-[di(ethyl-d₅)amino]-N-(2,6-dimethylphenyl)-

acetamide

MF: $C_{14}H_{12}D_{10}N_2O$

FW: 244.4

Chemical Purity: ≥98% (Lidocaine)


Deuterium

Incorporation: \geq 99% deuterated forms (d₁-d₁₀); \leq 1% d₀

Supplied as: A crystalline solid

Storage: -20°C Stability: ≥4 years

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.



Laboratory Procedures

Lidocaine-d₁₀ is intended for use as an internal standard for the quantification of lidocaine by GC- or LC-MS. The accuracy of the sample weight in this vial is between 5% over and 2% under the amount shown on the vial. If better precision is required, the deuterated standard should be quantitated against a more precisely weighed unlabeled standard by constructing a standard curve of peak intensity ratios (deuterated versus unlabeled).

Lidocaine- d_{10} is supplied as a crystalline solid. A stock solution may be made by dissolving the lidocaine- d_{10} in the solvent of choice, which should be purged with an inert gas. Lidocaine-d₁₀ is soluble in DMSO.

Description

Lidocaine is an inhibitor of voltage-gated sodium channels (Na_vs) and a local anesthetic.¹⁻³ It inhibits Na_v1.2 (K_i = 11 μ M), adult and neonatal Na_v1.5 (IC₅₀s = 380.1 and 360 μ M, respectively), and Na_v1.7 and Na_v1.8 channels (IC₅₀s = 450 and 104 μ M, respectively) expressed in Xenopus oocytes. Topical administration of lidocaine (1-2%) reduces escape behavior in a rhesus monkey model of noxious electrical cutaneous pain.⁴ Formulations containing lidocaine have been used as local and regional anesthetics.

References

- 1. Ragsdale, D.S., McPhee, J.C., Scheuer, T., et al. Common molecular determinants of local anesthetic, antiarrhythmic, and anticonvulsant block of voltage-gated Na+ channels. Proc. Natl. Acad. Sci. USA 93(17), 9270-9275 (1996).
- 2. Fraser, S.P., Onkal, R., Theys, M., et al. Neonatal Na, 1.5 channels: Pharmacological distinctiveness of a cancer-related voltage-gated sodium channel splice variant. Br. J. Pharmacol. (2021).
- Chevrier, P., Vijayaragavan, K., and Chahine, M. Differential modulation of Na, 1.7 and Na, 1.8 peripheral nerve sodium channels by the local anesthetic lidocaine. Br. J. Pharmacol. 142(3), 576-584 (2004).
- Lineberry, C.G. and Kulics, A.T. The effects of diazepam, morphine and lidocaine on nociception in rhesus monkeys: A signal detection analysis. J. Pharmacol. Exp. Ther. 205(2), 302-310 (1978).

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website

Copyright Cayman Chemical Company, 11/22/2022

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897

[734] 971-3335

FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.**CAYMANCHEM**.COM