



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Recombinant Streptavidin-NC

CLPRO338
CLPRO338-2
CLPRO338-3

Introduction: Streptavidin is a tetrameric protein secreted by *Streptomyces avidinii* which binds firmly to biotin. Streptavidin is widely used in molecular biology through its unique high affinity for the vitamin biotin. The dissociation constant (K_d) of the biotin-streptavidin complex is about ~10⁻¹⁵ mol/L. The strong affinity recognition of biotin and biotinylated molecules has made streptavidin one of the most important components in diagnostics and laboratory kits. The streptavidin/biotin system has one of the biggest free energies of association of yet observed for noncovalent binding of a protein and small ligand in aqueous solution (K_{assoc} = 10¹⁴). The complexes are also extremely stable over a wide range of temperature and pH.

Description: Recombinant Streptavidin-NC produced in *E. Coli* is a single, non-glycosylated, polypeptide chain having a molecular mass of 24kDa.

Source: *Escherichia Coli*.

Presentation: 20µg (CLPRO338), 100 µg (CLPRO338-2) or 1 mg (CLPRO338-2) sterile liquid formulation. The sterile solution contains 10mM K₂HPO₄-KH₂PO₄, PH 7.3.

Stability: rStreptavidin-NC although stable at 4°C for 3 weeks, should be stored below -18°C. Please prevent freeze thaw cycles.

Purity: Greater than 94.0% as determined by SDS-PAGE.

FOR RESEARCH USE ONLY

MW 06/05/20

Visit our website for your local distributor.

CEDARLANE®
www.cedarlanelabs.com

An ISO 9001:2015 and ISO 13485:2016
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com