



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7



[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 


# Falcon® Cell Scrapers

Improved designs to both the handle and the blade give Falcon Cell Scrapers superior handling and scraping ability



- Cross-ribbed polystyrene handle provides greater rigidity to ensure better control while scraping cells
- Highly compliant thermoplastic elastomer (TPE) blade easily pivots to provide multiple angles to remove cells from the entire growth surface
- Flexibility of the joint between the blade and handle improves ease of access into the neck of a flask or roller bottle
- Cell Scrapers are available in four sizes – choose the blade and handle lengths that best accommodate your needs
- Individually wrapped in peel-open, medical-style packaging for sterile presentation


The Cat. No. 353089 cell scraper has a larger blade to increase cell scraping efficiency in both Falcon dishes and 75 cm<sup>2</sup> flasks.



**FALCON®**

A Corning Brand

# Falcon® Cell Scrapers



- Larger handle diameter for better gripping
- Sharper scraping edge
- More compliant blade material


## Ordering Information

Cat. No.	Description	Qty.	Qty./Case
<b>Falcon Cell Scrapers</b>			
<b>1.8 cm highly compliant TPE blade</b>			
353085	18 cm polystyrene handle (recommended for use with Falcon 25 cm <sup>2</sup> flasks)	1	100
353086	25 cm polystyrene handle (recommended for use with Falcon 75 cm <sup>2</sup> flasks)	1	100
<b>3.0 cm highly compliant TPE blade</b>			
353087	40 cm polystyrene handle (recommended for use with large Falcon flasks)	1	100
353089	25 cm polystyrene handle (recommended for use with Falcon dishes and 75 cm <sup>2</sup> flasks)	1	100

To place an order in the U.S., contact Customer Service at:  
tel: **800.492.1110**, fax: 978.442.2476; email: **CLSCustServ@corning.com**.

For technical assistance, contact Technical Support at:  
tel: **800.492.1110**, fax: 978.442.2476; email: **CLSTechServ@corning.com**.

Outside the U.S., contact your local distributor or visit [www.corning.com/lifesciences](http://www.corning.com/lifesciences) to locate your nearest Corning office.



Unique blade design transfers scraping force to the blade edge, allowing it to conform to the scraping surface.



"Barbs" ensure blade remains attached to handle.

Corning acquired the Falcon® brand.  
For information, visit [www.corning.com/discoverylabware](http://www.corning.com/discoverylabware).

For a listing of trademarks, visit us at [www.corning.com/lifesciences/trademarks](http://www.corning.com/lifesciences/trademarks).  
Corning Incorporated, One Riverfront Plaza, Corning, NY 14831-0001

# CORNING

**Corning Incorporated**  
Life Sciences

836 North St.  
Building 300, Suite 3401  
Tewksbury, MA 01876  
t 800.492.1110  
t 978.442.2200  
f 978.442.2476

[www.corning.com/lifesciences](http://www.corning.com/lifesciences)