

Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

Elabscience®

Recombinant Mouse Siglec-E (C-Fc)

Catalog No. PKSM041446

Description	
Synonyms	SiglecE; Siglec-E; Sialic Acid Binding Ig-like Lectin E
Species	Mouse
Expression_host	Human Cells
Sequence	Gln20-Phe355
Accession	Q6PJ50
Mol_Mass	64.4 kDa
AP_Mol_Mass	85-95 kDa
Tag	C-Fc
Properties	
Purity	>95% as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg as determined by the LAL method.
Storage	Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 μ m filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.5.
Reconstitution	Please refer to the printed manual for detailed information.

Background

Siglecs are sialic acid specific I-type lectins that are characterized by an extracellular domain (ECD) with an N-terminal Ig-like V-set domain followed by varying numbers of Ig-like C2-set domains. Mouse Siglec-E, also known as Myeloid Inhibitory Siglec (MIS), is an 80 - 85 kDa member of the CD33-related subfamily of Siglecs. Rodent and primate Siglec gene families have significantly diverged, and Siglec-9 is the most likely human ortholog of mouse Siglec-E. Siglec-E is expressed as a heavily N-glycosylated disulfide-linked homodimer and shows binding preference for disialic acids in the alpha 2-8 linkage. Siglec-E is up-regulated and additionally phosphorylated following cellular stimulation by a variety of TLR agonists. Siglec-E signaling negatively regulates the LPS-induced production of TNF- alpha and IL-6 by macrophages. Its up-regulation in macrophages parallels the development of endotoxin tolerance. Siglec-E recognition of sialylated determinants on virulent T. cruzi contributes to the suppression of dendritic cell IL-12 p40 production.

Elabscience®

Tel:240-252-7368(USA) Fax:240-252-7376(USA) www.elabscience.com E-mail:techsupport@elabscience.com Elabscience Biotechnology Inc.

SDS-PAGE

