



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# SARS-CoV-2 (COVID-19) Spike FFPE 293T cell pellet slide

**Cat No. GTX435744**

<b>Application</b>	IHC-P (cell pellet)
<b>Species</b>	SARS Coronavirus 2

**Package**  
5 slide

## PRODUCT

### Summary

Cell pellet sections are 4 µm in thickness, and mounted on positively charged glass slides. Cell pellet consists of 293T cells overexpressing SARS-CoV-2 spike protein (full-length, no tag) fixed with 4% PFA (15 min, room temp) 48 hrs post-transfection, and then embedded in paraffin. Please see the corresponding overexpression cell lysate ([GTX535664](#)) and FFPE cell pellet block ([GTX435640](#)).

## APPLICATION

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Dilution
IHC-P (cell pellet)	Assay dependent


Not tested in other applications.

## PROPERTIES

**Storage** Store at room temperature.

**Note** For *In vitro* laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

## DATA IMAGES



### GTX435744 IHC-P (cell pellet) Image

SARS-CoV-2 (COVID-19) spike antibody detects SARS-CoV-2 (COVID-19) spike (full length) protein by immunohistochemical analysis.

Sample: Paraffin-embedded SARS-CoV-2 (COVID-19) Spike FFPE 293T cell pellet slide (GTX435744).

Green: SARS-CoV-2 (COVID-19) spike (full length) protein stained by SARS-CoV-2 (COVID-19) spike antibody (GTX135356) diluted at 1:1000.

Red: SARS-CoV-2 (COVID-19) spike (full length) protein stained by SARS-CoV / SARS-CoV-2 (COVID-19) spike antibody [1A9] (GTX632604) diluted at 1:1000.

Blue: Fluoroshield with DAPI (GTX30920).

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](#).