

Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

LW6

Cat. No.:	HY-13671		
CAS No.:	934593-90-	5	
Molecular Formula:	C ₂₆ H ₂₉ NO ₅		
Molecular Weight:	435.51		
Target:	HIF/HIF Prolyl-Hydroxylase; Apoptosis		
Pathway:	Metabolic Enzyme/Protease; Apoptosis		
Storage:	Powder	-20°C	3 years
		4°C	2 years
	In solvent	-80°C	6 months
		-20°C	1 month

®

MedChemExpress

SOLVENT & SOLUBILITY

DMF : 17.24 mg/mL	DMSO : 25 mg/mL (57.40 mM; Need ultrasonic) DMF : 17.24 mg/mL (39.59 mM; Need ultrasonic)				
		Solvent Mass Concentration	1 mg	5 mg	10 mg
	Preparing Stock Solutions	1 mM	2.2962 mL	11.4808 mL	22.9616 mL
		5 mM	0.4592 mL	2.2962 mL	4.5923 mL
		10 mM	0.2296 mL	1.1481 mL	2.2962 mL
	Please refer to the sol	lubility information to select the ap	opropriate solvent.		
In Vivo	1. Add each solvent one by one: 10% DMSO >> 40% PEG300 >> 5% Tween-80 >> 45% saline Solubility: 2.5 mg/mL (5.74 mM); Suspended solution; Need ultrasonic				
	2. Add each solvent one by one: 10% DMSO >> 90% corn oil Solubility: ≥ 2.5 mg/mL (5.74 mM); Clear solution				

BIOLOGICAL ACTIVITY		
Description	LW6 (HIF-1α inhibitor) is a novel HIF-1 inhibitor with an IC ₅₀ of 4.4 μM. LW6 decreases HIF-1α protein expression without affecting HIF-1β expression.	
IC ₅₀ & Target	IC50: 4.4 μM (HIF-1) ^[1]	
In Vitro	LW6 affects the stability of the HIF-1α protein. LW6 promotes the degradation of wild type HIF-1α, but not of a DM-HIF-1α with modifications of P402A and P564A, at hydroxylation sites in the oxygen-dependent degradation domain. LW6 induces the expression of von Hippel-Lindau (VHL), which interacts with prolyl-hydroxylated HIF-1α for proteasomal degradation. In the presence of LW6, knockdown of VHL does not abolish HIF-1α protein accumulation, indicating that LW6 degraded HIF-1α	

Product Data Sheet

OH

0

	via regulation of VHL expression ^[2] . In MDCKII-BCRP cells overexpressing BCRP, LW6 enhances significantly the cellular accumulation of mitoxantrone, a BCRP substrate. LW6 also down-regulates BCRP expression at concentrations of 0.1-10 μM ^[3] . LW6 inhibits the expression of HIF 1α induced by hypoxia in A549 cells at 20 μM, independently of the von Hippel Lindau protein. LW6 induces hypoxia selective apoptosis together with a reduction in the mitochondrial membrane potential ^[4] . MCE has not independently confirmed the accuracy of these methods. They are for reference only.
In Vivo	In mice carrying xenografts of human colon cancer HCT116 cells, LW6 demonstrates strong anti-tumor efficacy in vivo and causes a decrease in HIF-1α expression in frozen-tissue immunohistochemical staining ^[2] . MCE has not independently confirmed the accuracy of these methods. They are for reference only.

PROTOCOL

Cell Assay ^[2]	Inhibition of HIF-1a is assayed by a reporter assay using dualluciferase reporter assay system. HCT116 cells in 75-90% confluence are transiently co-transfected with pGL3-HRE-luciferase plasmid containing six copies of HREs from human VEGF genes and pRLSV40 encoding firefly renilla luciferase and incubated for 24 h. Cells are treated with LW6 or 17-AAG for 16 h before report assay. Luciferase activity is integrated over a 10 second period and measured using a luminometer ^[2] . MCE has not independently confirmed the accuracy of these methods. They are for reference only.
Animal Administration ^[2]	Mice: The mice receive the following treatments using a dosing vehicle solution, containing 10% dimethylacetamide, 10% Cremophor EL and 80% of sodium carbonate buffer (pH 10), by intraperitoneal (i.p.) injection: group1(control group; six mice), vehicle solution; group2 (six mice), LW6 at a dose of 10 and 20mg/kg (QD); and group 3 (six mice), topotecan at a dose of 2mg/kg, (Q2D), which is the dose and dosing schedule that showed more than 60% inhibition of growth of HCT116 tumors. The treatments are continued for 13 days ^[2] . MCE has not independently confirmed the accuracy of these methods. They are for reference only.

CUSTOMER VALIDATION

- Cell Metab. 2023 Jun 14;S1550-4131(23)00209-7.
- Cell Death Dis. 2021 May 14;12(5):490.
- Oncogene. 2022 Nov 18.
- iScience. 2021 May 7;24(6):102521.
- Int J Mol Sci. 2023 Jun 16, 24(12), 10236.

See more customer validations on <u>www.MedChemExpress.com</u>

REFERENCES

[1]. Naik R, et al. Synthesis and structure-activity relationship study of chemical probes as hypoxia induced factor-1α/malate dehydrogenase 2 inhibitors. J Med Chem. 2014 Nov 26;57(22):9522-38.

[2]. Lee K, et al. LW6, a novel HIF-1 inhibitor, promotes proteasomal degradation of HIF-1alpha via upregulation of VHL in a colon cancer cell line. Biochem Pharmacol. 2010 Oct 1;80(7):982-9.

[3]. Song JG, et al. Discovery of LW6 as a new potent inhibitor of breast cancer resistance protein.

[4]. Sato M, et al. LW6, a hypoxia-inducible factor 1 inhibitor, selectively induces apoptosis in hypoxic cells through depolarization of mitochondria in A549 human lung cancer cells. Mol Med Rep. 2015 Sep;12(3):3462-8.

Caution: Product has not been fully validated for medical applications. For research use only.

 Tel: 609-228-6898
 Fax: 609-228-5909
 E-mail: tech@MedChemExpress.com

 Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA